# GIARDIA CEL

English  
Product Code: KG1

FIGURE 1: GIARDIA CEL DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE GIARDIA CEL

TABLEAU 1: SENSIBILITÉ, SPÉCIFICITÉ ET AUTRES DONNÉES DU TEST GIARDIA CEL  
 TABELLE 1: SENSITIVITÄT, SPEZIFIZITÄT UND ANDERE DATEN ZUM GIARDIA CEL  
 TABELLA 1: SENSIBILITA', SPECIFICITA' ED ALTRI DATI SULLA GIARDIA CEL  
 TABLA 1: SENSIBILIDAD, ESPECIFICIDAD Y OTROS DATOS DEL GIARDIA CEL  
 TABELA 1: SENSIBILIDADE, ESPECIFICIDADE, & OUTROS DADOS DO GIARDIA CEL

Trial Essai Versuch Prova Prueba Teste	Sensitivity Sensibilité Sensitivität Sensibilita' Sensibilidad Sensibilidade	Specificity Spécificité Spezifität Specificita' Especificidad Especificidade	Repeatability Répétabilité Wiederholpräzision Ripetibilità Repetibilidad Repetição	Reproducibility Reproductibilité Reproduzierbarkeit Riproducibilità Reproducibilidad Reprodutibilidade
A	100%	100%	-	-
B	94%	100%	-	-
C	100%	100%	-	-
D	100%	100%	-	-
E	90.3%	100%	-	-
F	-	-	100% Correlation	100% Correlation

**Not cross reactive with /** Pas de Réaction Croisée avec / Keine Kreuzreaktionen mit / Non mostra reazione crociata con / No muestra reacción cruzada con / Não apresenta reações cruzadas com:  
*Strongyloides stercoralis* larva, *Enterobius vermicularis* ova, *Opisthorchis viverrini* ova, Hook Worm ova, *Trichuris trichuria* ova, *Hymenolepis nana* ova, *Blastocystis hominis* cysts, *Entamoeba coli* cysts, *Endolimax nana* cysts, *Entamoeba hartmanni* cysts, *Trichomonas hominis* cysts

## EXPLANATION OF SYMBOLS

- Consult Instructions for Use
- In Vitro Diagnostic Medical Device
- Temperature Limitation
- Batch
- Control Positive
- Use By/Expiration Date
- Do Not Re-use

Cellabs Pty Ltd  
 Unit 7, 27 Dale Street  
 Brookvale, NSW 2100 Australia  
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426  
 Web: <http://www.cellabs.com.au>  
 Email: [sales@cellabs.com.au](mailto:sales@cellabs.com.au)

WMDE  
 Bergerweg 18  
 6085 AT Horn  
 The Netherlands

LG1.12.1  
 23 May 2017

## INTENDED USE AND PRINCIPLE OF THE TEST

The Giardia Cel IF Test is an *in vitro* direct immunofluorescence test for the detection of the protozoan parasite *Giardia* in faecal and environmental specimens. The fluorescein-labelled mouse monoclonal antibody reagent binds specifically to *Giardia* cysts in the specimen. *Giardia* cysts display bright green fluorescence with typical morphology.

## CONTENTS OF THE KIT

		KG1	Bulk
	Giardia Cel Reagent	1.25mL	5mL
	Positive Control Slide (Single Use only)	1	-
	Mounting Fluid	2.5mL	-
	Tests	50	200

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

## MATERIALS REQUIRED BUT NOT PROVIDED

Microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; acetone for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x200-x400 magnification.

## PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix components from different kits. The Giardia Cel Reagent has been optimised for use with Cellabs Positive Control Slide and Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. The detection of *Giardia* from water is dependent on the efficient filtration and recovery from large volumes of water (500-1000L). A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

## SPECIMEN PREPARATION

Stool specimens, either fresh or preserved in 10% formalin or SAF may be used. A procedure for processing fresh stool specimens is described below. If stool specimens are not processed immediately they may be stored at 2-8°C for 24-48 hours. Stools may be concentrated prior to processing using concentration kits. They may also be concentrated by the formalin-ether sedimentation method.

Prepare an approximate 1/10-1/50 dilution of stool. Dilute 50µL or 5mm diameter of stool in PBS with 0.1% sodium azide and mix thoroughly to disperse the specimen using an applicator stick and allow large stool particles to sediment.

## PREPARATION OF SLIDES

- Place 20µL of the faecal or concentrated specimen from water onto a microscope slide.
- Allow the specimen to completely air dry.
- Fix the slides for 5 minutes in acetone and allow to air dry.

## INSTRUCTIONS FOR USE

- Add 25µL of to the fixed specimen and , covering well area.
- Incubate the slides at 37°C in a humid chamber for 30 minutes. Do not allow the slides to dry as this may cause non-specific binding.
- Rinse gently in a bath of PBS for one minute.
- Drain slide and remove excess moisture around well with tissue.
- Add a drop of to the slide well. Place a coverslip on top of the drop and remove air bubbles.
- Scan the entire specimen using a fluorescence microscope initially at X200 magnification, then at X400 for confirmation. Read immediately or store at 2-8°C in the dark for up to 24 hours.

## READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

*Giardia* cysts are 8-12µm in size, elliptical in shape with bright green fluorescence. A test may be considered positive if one or more cysts are present. In active infections the cyst numbers are moderate to high, while in the asymptomatic stage, cyst numbers are generally low. If the test is negative and symptoms continue, another specimen should be considered as cyst excretion may vary from day to day. A comparison should be made with the positive control slide for the appearance of organisms. Negative specimens should have signs of faecal material containing debris such as undigested food particles, cells (counterstained red-brown), and possibly other organisms such as bacteria and protozoa, and should not show any organisms of characteristic green fluorescence and morphology. Non-specific fluorescence may occur in some specimens; however this may be distinguished morphologically from the *Giardia* cysts. Algae may be present in environmental specimens and these fluoresce red-orange.

## WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

## SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE GIARDIA CEL

Refer to summary table at end of insert. All data on the Giardia Cel can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

## INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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