GIARDIA CELISA

INTENDED USE AND PRINCIPLE OF THE TEST

The Giardia CELISA kit is a qualitative *in vitro* enzyme immunoassay for the detection of *Giardia lamblia* cyst antigens in faecal specimens. The Giardia CELISA detects *Giardia* antigens in a capture enzyme immunoassay. Antigens from faecal specimens are bound to microplates, which have been coated with purified mouse monoclonal antibody to *Giardia*. After a washing step, the conjugate (consisting of rabbit polyclonal antibody) is added, and binds to the antibody/antigen complex. After another washing step, the chromogenic substrate 3, 3', 5, 5'-tetramethyl benzidine (TMB) is added, causing a colour change. Enzyme activity is directly proportional to the concentration of *Giardia* antigen present in the sample and control.

CONTENTS OF THE KIT

GMW	Celisa Plate – 1x96 wells - (single use only)	1 plate
CONTROL +	Positive Control (Black Cap)	3.5mL
GDB	Diluent (White Cap)	40mL
GPO	Enzyme Conjugate (Red Cap)	7.0mL
GWB	Wash Buffer Concentrate (20x) (White Cap)	50mL
GSC	Substrate Solution (Blue Cap)	14.0mL
GSS	Stopping Solution (Yellow Cap)	7.0mL
	Graduated Disposable Pipettes	100
	Plastic Adhesive Sheets	1
	Resealable Bag	1

Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Distilled water for diluting wash reagent; 1L bottle for diluting wash reagent; Squirt bottle for wash reagent; Test tubes; Applicator sticks; Micropipettes with disposable tips; Vortex mixer; Timer; Discard container and absorbent paper; and Spectrophotometer capable of reading absorbance of EIA plates at 450nm or 450/620nm.

PRECAUTIONS

For *in vitro* diagnostic use only. Reagents should not be used after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Caps and tips are colour coded, do not mix. When handling assay wells, avoid scratching the bottom of the wells as this may result in elevated absorbance readings. Thimerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eyes and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. The positive control has been shown to be non-infectious in tissue culture; however it should be handled as though potentially infectious. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Material Safety Data Sheet.

INSTRUCTIONS FOR USE

Preparation of Wash Buffer

If crystals are present in the concentrate, warm to dissolve. Prepare working strength **WASH BUFFER** by adding 50mL of **GWB** to 950mL of distilled water. Label the bottle and store at 2-8°C.

Specimen Collection and Preparation

Fresh/Frozen Specimens: Frozen specimens should be thawed. Prepare **SAMPLE** by adding 0.1mL specimen to 0.4mL **GDB** to make a 1/5 dilution. If the specimen cannot be pipetted, use an applicator stick to transfer approximately 0.1g of specimen (about the size of a small pea 4mm in diameter).

Preserved Specimens: To prepare **SAMPLE** mix contents of container thoroughly before transferring specimen. No further processing or dilution is required.

Assay Procedure

- Bring all reagents to room temperature (18-25 °C) before use.
- Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
- Remove required number of **GMW** strips. Reseal the foil bag containing unused microwell strips immediately with tape.
- Two (2) control wells must be used each time the test is performed. These serve as positive and negative controls. Add one (1) drop (50µL) of **CONTROL +** to the positive control well and two (2) drops (100µL) of Diluent **GDB** to the negative control well.
- Add two (2) drops (100µL) of **GDB** to the test well(s). Using plastic pipettes, transfer (1) drop (50µL) of **SAMPLE** to well(s) containing **GDB**. Tap to mix, seal control and test wells and incubate at room temperature for one (1) hour.
- Shake out the contents of the assay wells. Wash the wells by using the **WASH BUFFER** in a squirt bottle with fine-tipped nozzle, directing the wash buffer to the bottom of the well with force. Fill the wells, and then shake the wash buffer out. Bang the inverted plate on a dry paper towel and repeat the washing step another three (3) times. If any particulate matter is in the wells continue to wash until it is removed.
- After washing, bang dry (until no liquid comes out) and dispose of paper towels and specimen containers.
- Add one (1) drop (50µL) of Enzyme Conjugate **GPO** (red cap) to each well, tap to mix, seal control and test wells, and incubate for 30 minutes at room temperature.
- Repeat the washing and steps as before (Steps 6 and 7).
- Add two (2) drops (100µL) of Substrate Solution **GSC** (blue cap) to each well. Gently tap the wells to mix the substrate. Incubate for ten (10) minutes at room temperature in the dark. Tap the wells two (2) or three (3) times during incubation to mix the contents. Observe the blue colour formation.
- Dispense one (1) drop (50µL) of Stopping Solution **GSS** (yellow cap) to each well. The addition of the stop solution changes the colour of the wells from blue to yellow.
- Quantitate the results by measuring the absorbance in a spectrophotometer at a wavelength of 450nm (reference 620nm on dual wavelength reader). The instrument should be blanked against air. Wipe the underside of each well before measuring the optical density. Plates should be read within 10 minutes of stopping the reaction. If a spectrophotometer is unavailable, the test may be read visually in good light against a white background (such as a sheet of paper).

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

Visual

If *Giardia* antigens are present in the stool specimens, the positive wells will show a yellow colour that is obviously more yellow than the negative control well. Negative wells will have the same or less colour than the negative control well.

Spectrophotometer Readings

	Negative Control	Positive Control	Negative Specimen	Positive Specimen
Single Wavelength at 450nm	< 0.150	≥0.500	< 0.150	≥0.150
Dual Wavelength at 450/620nm	<0.090	≥0.500	<0.090	≥0.090

The Positive and Negative Control must fall in their respective ranges or the test is not valid. A positive test result indicates that *Giardia* cyst antigen is present in the faecal specimen. A negative result indicates that *Giardia* cyst antigen is not present in the faecal specimen.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA

Refer to summary table at end of insert. All data on the Giardia CELISA can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances

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Cellabs Pty Ltd