TOXOCARA IgG CELISA

INTENDED USE AND PRINCIPLE OF THE TEST

This test is designed to detect and measure antibodies to *Toxocara canis* excretory-secretory (ES) antigen produced during infection in humans. The Toxocara IgG CELISA is a solid phase enzyme immunoassay. Antibody is bound by the *Toxocara* ES antigen coated on the inner surface of the test strip. Peroxidase-conjugated antibody to human IgG is added and reacts with bound antibody. A chromogenic substrate for peroxidase is added. If antibody to *Toxocara* is present there is a reaction which results in the development of a blue colour, which is in proportion to the serum level of the antibody. Addition of stopping solution changes the blue colour to yellow.

CONTENTS OF THE KIT

TXMW	Celisa Plate – 1x96 wells - (single use only)	2 plates
CONTROL +	Positive Control	0.2mL
CONTROL -	Negative Control	0.2mL
TXPO	Enzyme Conjugate (200x)	0.2mL
TXPT	PBS/Tween (20x)	110mL
TXSC	Substrate Chromogen (TMB) (20x)	1.5mL
TXSB	Substrate Buffer	24mL
TXSS	Stopping Solution	12mL

Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Micropipettes and tips; clean glassware or plastic containers for solutions; humid chamber; ELISA washer; and Spectrophotometer to read absorbances at a single wavelength of 450nm.

PRECAUTIONS

For in vitro diagnostic use only. Reagents should not be used after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Thimerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eyes and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Material Safety Data Sheet.

INSTRUCTIONS FOR USE

Preparation of Wash Buffer

If crystals are present, warm the concentrate to dissolve. For each microplate, add 50mL PBS-Tween concentrate **TXPT** to 950mL of distilled water. Label the bottle **WASH BUFFER**. Store at 2-8°C.

Preparation of Samples

Collect patient blood samples by standard venipuncture procedure. The test may be used with serum or plasma. The serum or plasma should be stored below -10°C if the analysis is delayed. Prepare a 1/100 dilution of the positive control **CONTROL +**, the negative control **CONTROL -** and the test (patient specimen) sera samples in **WASH BUFFER**, ensuring proper mixing. Record the position of each diluted **SAMPLE** on a work sheet.

Assay Procedure

- Bring all reagents to room temperature (18-25 °C) before use.
- Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
- Remove required number of **TXMW** strips. Reseal the foil bag containing unused microwell strips immediately with tape.
- Pipette 100µL of the **SAMPLE**, **CONTROL +** and **CONTROL -**, into individual microwells. Include two positive and two negatives in each assay run. Cover and incubate for one (1) hour at room temperature (RT) in a humid chamber.
- In the last 10 minutes of the incubation period, prepare the working strength **CONJUGATE**. Add 5µL of Enzyme Conjugate **TXPO** to 995µL of **WASH BUFFER** and mix thoroughly (allow 1mL per strip of 8 wells).
- Wash the wells preferably using an automatic plate/strip washer or manually as follows:
 - Empty contents from the wells. Refill with the **WASH BUFFER**.
 - Repeat this process a further three (3) times. After the fourth wash, bang inverted wells dry on absorbent tissue.
 - NB: take care when flicking out plates, hold side of frame firmly to hold strips in place.
- Add 100µL of **CONJUGATE** to each well. Incubate for one (1) hour at RT in a humid chamber.
- In the last 10 minutes of the incubation period, prepare the working strength **SUBSTRATE**. Add 50µL of Substrate Chromogen **TXSC** to 950µL of Substrate Buffer **TXSB** and mix thoroughly (allow 1mL per strip of 8 wells). The stability of the solution is 30 minutes.
- Repeat washing as in step 6.
- Add 100µL of fresh **SUBSTRATE** and incubate in the dark (covered) at room temperature for 15 minutes.
- Add 50µL of Stopping Solution **TXSS**. Tap the plate to mix.
- Read the results visually or in a spectrophotometer at 450nm.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

Photometrically

It is optimal to use a dual beam ELISA reader. Blank machine against the well designated as the blank. Within 30 minutes of stopping the reaction read absorbance of all wells including all the controls at 450nm on an ELISA reader.

For the test results to be accepted the controls must read as follows:

Control	O.D Value (450nm)
Negative Control	< 0.25
Positive Control	> 1.0

Test results should be interpreted as follows:

O.D Value (450nm)	Interpretation
< 0.25	NEGATIVE specimen
0.25 - < 0.5	Contains antibodies - However, level is lower than generally accepted significant level
0.5 - 1.0	Significant antibody level - Indicative of past infection or a current low level infection
> 1.0	Strong POSITIVE result - Indicative of recent clinical Toxocariasis

Visual Reading

Results can also be read visually. Note colour of the specimen wells and compare with the colour of the negative control. If the colour, read visually, of a specimen well is above that of the negative control, then the specimen contains antibody.

A positive result usually indicates that the donor is or has been infected with *Toxocara canis*. It does not necessarily mean that the donor has clinical Toxocariasis. The test may be negative early in the infection.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE TOXOCARA CELISA

Refer to summary table at end of insert. All data on the Toxocara CELISA can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances

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FIGURE 1 TOXOCARA CELISA DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE TOXOCARA IgG CELISA
 TABLEAU 1: SENSIBILITE, SPECIFICITE ET AUTRES DONNEES DU TEST TOXOCARA IgG CELISA
 TABELLE 1: SENSITIVITÄT, SPEZIFITÄT UND ANDERE DATEN ZUM TOXOCARA IgG CELISA
 TABLA 1: SENSIBILIDAD, ESPECIFICIDAD Y OTROS DATOS DEL TOXOCARA IgG CELISA
 TABELA 1: SENSIBILIDADE, ESPECIFICIDADE E OUTROS DADOS DO TOXOCARA IgG CELISA

Trial Essai Versuch Prueba Teste	Sensitivity Sensibilité Sensitivität Sensibilidad Sensibilidade	Specificity Spécificité Spezifität Especificidad Especificidade	Repeatability Répétabilité Wiederholpräzision Repetibilitàd Repetição	Reproducibility Reproductibilité Reproduzierbarkeit Reproducibilidad Reproduzibilidade
A			Positive CV = 7.11% Negative CV = 9.36%	Positive CV = 7.24% Negative CV = 9.49%
B	90%	94%		

EXPLANATION OF SYMBOLS

	Consult Instructions for Use		Cellabs Pty Ltd Unit 7, 27 Dale Street Brookvale, NSW 2100 Australia Tel: +61 2 9905 0133 Fax: +61 2 9905 6426 Web: http://www.cellabs.com.au Email: sales@cellabs.com.au
	In Vitro Diagnostic Medical Device		WMDE Bergeweg 18 6085 AT Horn The Netherlands
	Temperature Limitation 2°C - 8°C		Insert Version LT3.12.1 23 May 2017
	Batch		
	Control Positive		
	Control Negative		
	Use By/Expiration Date		
	Do Not Re-use		