# MALARIA Ag CELISA

## INTENDED USE AND PRINCIPLE OF THE TEST

The Malaria Ag CELISA™ kit has been designed as a confirmatory test for *falciparum* malaria in situations where traditional diagnosis is unclear, for screening blood transfusion products, or to confirm cases of travel-related infection. It is **not** intended to replace the conventional blood film diagnosis. The sandwich ELISA principle is employed. The microwells are pre-coated with anti-*P. falciparum* monoclonal capture antibody. A conjugate of enzyme labelled anti-*P. falciparum* monoclonal antibody is prepared and incorporated into the kit. The user adds a blood sample to the coated wells and if *P. falciparum* malaria antigen is present it binds to the coated well. All other blood components are removed by a washing step. The horse radish peroxidase enzyme labelled anti-malaria monoclonal indicator antibody conjugate is then added. It binds to any *falciparum* malaria antigen that has been captured on the well surface. The strip is then washed and the enzyme substrate solution is then added to the test wells and incubated. Colour generated indicates that *P. falciparum* malarial antigen in the test sample.

## CONTENTS OF THE KIT

The Malaria Ag CELISA is available in 3 different formats:

		Standard	Automated	Bulk
<b>MAMW</b>	Celisa Plate – 1x96 wells - (single use only)	2 plates	2 plates	10 plates
<b>MAPC</b>	Positive Control	-	1 x 0.5mL	-
<b>MANC</b>	Negative Control	-	1 x 0.5mL	-
<b>MAPO</b>	Enzyme Conjugate (200x)	1 x 0.12mL	1 x 0.24mL	1 x 1.2mL
<b>MACD</b>	Conjugate Diluent	1 x 24mL	2 x 18mL	1 x 120mL
<b>MAPT</b>	PBS/Tween (20x)	1 x 125mL	1 x 125mL	1 x 650mL
<b>MASC</b>	Substrate Chromogen (TMB) (20x)	1 x 1.2mL	1 x 2.4mL	1 x 7.5mL
<b>MASB</b>	Substrate Buffer	1 x 24mL	2 x 18mL	1 x 125mL
<b>MASS</b>	Stopping Solution	1 x 12mL	1 x 18mL	2 x 30mL

Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

## MATERIALS REQUIRED BUT NOT PROVIDED

Malaria positive blood samples, distilled water, micropipettes and tips, clean glassware or plastic containers for solutions, humid chamber, ELISA washer, spectrophotometer to read absorbances at a single wavelength of 450nm, or at dual wavelengths of 450nm and 620nm.

## PRECAUTIONS

For in vitro diagnostic use only. Reagents should not be used after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Thimerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eyes and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Material Safety Data Sheet.

## INSTRUCTIONS FOR USE

### Preparation of Wash Buffer

If crystals are present, warm the concentrate to dissolve. For each microplate, add 50mL PBS-Tween concentrate **MAPT** to 950mL of distilled water. Label the bottle **WASH BUFFER**. Store at 2-8°C.

### Preparation of Samples

Collect patient blood by standard venipuncture procedure using an anticoagulant. Lyse blood by freezing, use the lysed blood as the test specimen **SAMPLE**. Serum or plasma may be used as an alternative to whole blood but the use of these samples may result in the loss of sensitivity. The blood should be stored below -10°C if the analysis is delayed.

### Positive and Negative Controls

For Standard and Bulk versions of the kit, use the Conjugate Diluent **MACD** undiluted as the kit Negative Control, and use a confirmed Malaria positive blood sample as the Positive Control.

### Assay Procedure

1. Bring all reagents to room temperature (18-25 °C) before use.
2. Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
3. Remove required number of **MAMW** strips. Reseal the foil bag containing unused microwell strips immediately with tape.
4. Pipette 100µL of the **SAMPLE**, positive control (**MAPC** or known positive) and negative control (**MANC** or **MACD**), into individual microwells. Include two positive and two negatives in each assay run. Cover and incubate for one (1) hour at room temperature in a humid chamber.
5. In the last 10 minutes of the incubation period, prepare the working strength **CONJUGATE**. Add 5µL of conjugate concentrate **MAPO** to 995µL of conjugate diluent **MACD** and mix thoroughly (allow 1mL per strip of 8 wells).
6. Wash the wells preferably using an automatic plate/strip washer or manually as follows:
  - Empty contents from the wells. Refill with the **WASH BUFFER**.
  - Repeat this process a further four (4) times. After the fifth wash, bang inverted wells dry on absorbent tissue.
  - NB: take care when flicking out plates, hold side of frame firmly to hold strips in place.
7. Add 100µL of **CONJUGATE** to each well. Incubate for one (1) hour at room temperature in a humid chamber.
8. In the last 10 minutes of the incubation period, prepare the working strength **SUBSTRATE**. Add 50µL of substrate chromogen **MASC** to 950µL of substrate buffer **MASB** and mix thoroughly (allow 1mL per strip of 8 wells). The stability of the solution is 30 minutes.
9. Repeat washing as in step 6.
10. Add 100µL of fresh **SUBSTRATE** and incubate in the dark (covered) at room temperature for 15 minutes.
11. Add 50µL of Stopping Solution **MASS**. Tap the plate to mix.
12. Read the results visually or in a spectrophotometer at 450nm, or 450nm/620nm, blanking the machine on air.

## READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

### Visually

Observe the colour intensity of the control and specimen wells. The Positive Control should be blue before, and yellow after stopping.

### Photometrically

Read the microwell plate at 450nm or 450nm / 620nm in a compatible ELISA plate reader, blanked against air.

For the test results to be accepted the Negative Control must read as follows:

	O.D Value (450nm, 450/620nm)
Negative Control	< 0.1
Cut-Off level	= Negative Control OD + 0.1

Negative blood samples should give optical density readings below 0.1 OD units. However, to allow for inter-laboratory variation we strongly recommend that each laboratory run a number of known negative blood samples to allow standardisation of the CELISA positive / negative cut-off level.

All specimens with an absorbance value above the cut-off level should be considered positive for *P. falciparum* antigen. A positive result indicates the presence of *P. falciparum* antigen. This is suggestive of current or very recent infection. The assay has been shown to detect *P. falciparum* infection at parasitaemias as low as 0.001%. The intensity of colour is not proportional to the level of parasitaemia. *P. vivax*, *P. ovale* and *P. malariae* infections are not detected. Please note that the test may remain positive for several days after parasites are no longer detectable in blood films.

#### WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

#### SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE MALARIA AG CELISA

Refer to summary table at end of insert. All data on the Malaria Ag CELISA can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

#### INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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MDCI Ltd
Cellabs Pty Ltd