PNEUMO CEL

FIGURE 1: PNEUMO CEL DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE PNEUMO CEL

TABEAU 1: SENSIBILITÉ, SPÉCIFICITÉ ET AUTRES DONNÉES DU TEST PNEUMO CEL
 TABELLE 1: SENSITIVITÄT, SPEZIFITÄT UND ANDERE DATEN ZUM PNEUMO CEL
 TABELLA 1: SENSIBILITA', SPECIFICITA' ED ALTRI DATI SULLA PNEUMO CEL
 TABLA 1: SENSIBILIDAD, ESPECIFICIDAD Y OTROS DATOS DEL PNEUMO CEL
 TABELA 1: SENSIBILIDADE, ESPECIFICIDADE, & OUTROS DADOS DO PNEUMO CEL

| Trial Essai Versuch Prova Prueba Teste | Sensitivity Sensibilité Sensitivität Sensibilita' Sensibilidad Sensibilidade | Specificity Spécificité Spezifität Specificita' Especificidad Especificidade | Repeatability Répétabilité Wiederholpräzision Ripetibilita Repetibilidad Repetição | Reproducibility Reproductibilité Reproduzierbarkeit Riproducibilita Reproducibilidad Reprodutibilidade |
|---|---|---|---|---|
| A | 100% | 100% | - | - |
| B | 100% | 100% | - | - |
| C | 100% | 100% | - | - |
| D | 100% | 95.1% | - | - |
| E | - | - | 100% Correlation | 100% Correlation |
| Not cross reactive with / Pas de Réaction Croisée avec / Keine Kreuzreaktionen mit / Non mostra reazione crociata con / No muestra reacción cruzada con / Não apresenta reações cruzadas com: <i>Cryptococcus neoformans, Candida tropicalis, Saccharomyces cerevesiae, Aspergillus niger, Candida krusei, Candida albicans</i> | | | | |

EXPLANATION OF SYMBOLS

Consult Instructions for Use
 Consulter le Mode d'Emploi
 Gebrauchsanweisung
 Consultare le istruzioni per l'uso
 Consultar el manual de instrucciones
 Consultar o Folheto de Instruções

In Vitro Diagnostic Medical Device
 Produit Diagnostique Médical In Vitro
 In-vitro-Diagnostika
 Diagnostico in vitro Dispositivo medico
 Producto Sanitario para Diagnóstico In Vitro
 Produto de Diagnóstico Medico In vitro

Temperature Limitation
 Températures Limites
 Temperaturbegrenzung
 Temperatura limite
 Limite de temperatura
 Temperatura Limitada

Biological Risks
 Risques Biologiques
 Biogefährdung
 Rischi biologici
 Riesgo biológico
 Riscos Biológicos

Manufacturer
 Fabricant
 Hersteller
 Fabricante
 Fabricante
 Fabricante

Cellabs Pty Ltd
 Unit 7, 27 Dale Street (PO Box 421)
 Brookvale, NSW 2100 Australia
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426
 Web: <http://www.cellabs.com.au>
 Email: sales@cellabs.com.au

Authorised Representative in the EC
 Représentant Autorisé dans la EC
 Bevollmächtigter in der Europäischen Gemeinschaft
 Rappresentante Autorizzato nella EC
 Representante autorizado en la EC
 Representante Autorizado na CE

MDCI Ltd
 Arundel House, 1 Liverpool Gardens
 Worthing, West Sussex BN11 1SL
 United Kingdom

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LP1.11 - PDF
 4th December 2006

INTENDED USE AND PRINCIPLE OF THE TEST

The Pneumo Cel IF Test is a rapid *in vitro* direct immunofluorescence test for the detection of *Pneumocystis carinii* in human specimens. The fluorescein-labelled mouse monoclonal antibody reagent binds specifically to *P. carinii* in the specimen. *P. carinii* display bright green fluorescence with typical morphology.

CONTENTS OF THE KIT

| | RP1 | RP | RMF | KP1 Demo 0.25mL | KP1 Standard 1.25mL | Bulk 5mL |
|--|-----|----|-----|--------------------|------------------------|-------------|
| Pneumo Cel Reagent | | | | 1 | 1 | - |
| Positive Control Slide (Single Use only) | | | | 2.5mL | 2.5mL | - |
| Mounting Fluid | | | | 10 | 50 | 200 |
| Tests | | | | | | |

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; acetone for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x400-x1000 magnification.

PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. The Pneumo Cel Reagent has been optimised for use with Cellabs Positive Control Slide and Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

SPECIMEN COLLECTION

Clinical specimens should be processed as soon as possible after collection. Induced sputum, bronchial wash, bronchoalveolar lavage and tissue impression smears may be used. The presence of mucous in specimen may prevent adequate staining. Work in a biohazard hood. A positive control slide should be included with each series of test specimens.

Preparation of Induced Sputum Specimens

Prepare 0.3% dithiothreitol (DTT) in distilled water. Mix equal volumes of sputum and DTT solution. Mix vigorously (vortex) to dispense clumps. Incubate for 30 minutes at 37°C. Centrifuge for 5 minutes at 1500g. Decant supernatant. Resuspend pellet in 0.2-0.5mL of PBS and mix vigorously

Preparation of Bronchial Wash and Bronchoalveolar Lavage Specimens

Centrifuge the specimen at 1500g for 10 minutes. Remove the supernatant leaving 0.5mL. Decant supernatant. Resuspend pellet in 0.2-0.5mL of PBS and mix vigorously

PREPARATION OF SLIDES

Apply 25µL to 1 well (0.8cm) slide and air dry. Fix slide for 5 minutes in acetone. If slides are not stained immediately they may be stored below 0°C for several weeks.

INSTRUCTIONS FOR USE

- Add 25µL of RP1 to the fixed specimen and RP, covering well area.
- Incubate the slides at 37°C in a humid chamber for 30 minutes. Do not allow the slides to dry as this may cause non-specific binding.
- Rinse gently in a bath of PBS for one minute.
- Drain slide and remove excess moisture around well with tissue.
- Add a drop of RMF to the slide well. Place a coverslip on top of the drop and remove air bubbles.
- Scan the entire specimen using a fluorescence microscope initially at X400 magnification, then at X1000 for confirmation. Read immediately or store at 2-8°C in the dark for up to 24 hours.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

In contrast to the reddish-brown colour of counter-stained material, the pleomorphic trophozoites of *Pneumocystis carinii* (PC) exhibit an intense apple-green fluorescence, while the thick-walled cysts appear as round bodies stained at the rim in a "halo-like" appearance. The organisms occur either singularly or more commonly in clusters in an extracellular matrix, which is also stained. A test is considered positive if two or more cysts are present. The control specimen should be used as a comparison for the size and appearance of PC organisms. Negative specimens should not display fluorescent organisms of characteristic morphology. Fluorescent material which is not apple green and can be distinguished from typical PC forms should be disregarded.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE PNEUMO CEL

Refer to summary table at end of insert. All data on the Pneumo Cel can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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