CRYPTO/GIARDIA CEL

FIGURE 1: CRYPTO/GIARDIA CEL DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE CRYPTO/GIARDIA CEL

Trial Essai Versuch TABELLA 1: TABELLA 1: TABELLA 1: TABELLA 1: TABELLA 1:	Sensitivity Sensibilité Sensitivität Sensibilität Sensibilità Sensibilidat Sensibilidade Sensibilitet	Specificity Spécificité Spezifität Specificita' Especificidad Especificidade Specificitet	Repeatability Répétabilité Wiederholpräzision Wiederholbarkeit Ripetibilita Repetibilitad Repetição Repetierbarhed	Reproducibility Reproductibilité Reproduzierbarkeit Riproducibilita Reproducibilidat Reprodutibilidade Reproducerbarhed
A	100%	100%	-	-
B	-	-	100% Correlation	100% Correlation
<p>Not cross reactive with / Pas de Réaction Croisée avec / Keine Kreuzreaktionen mit / Non mostra reazioni crociata con / No muestra reacción cruzada con / Não apresenta reações cruzadas com / Krydsreagerer ikke med: <i>Strongyloides stercoralis larva, Enterobius vermicularis ova, Opisthorchis viverrini ova, Hook Worm ova, Trichuris trichiura ova, Hymenolepis nana ova, Blastocystis hominis cysts, Entamoeba coli cysts, Entamoeba hartmanni cysts, Eimeria tenella, Toxoplasma gondii, Trichomonas hominis, Iospora belli, Microsporidium sp., Pneumocystis carinii, Candida sp., Aspergillus sp., Listeria sp., Legionella sp., Erysipelothrix sp., Escherichia sp., Lactobacillus sp., Streptococcus sp., Staphylococcus sp., Bacillus sp., Clostridium sp., Brochothrix sp., Corynebacterium sp.</i></p>				

EXPLANATION OF SYMBOLS

Consult Instructions for Use
 Consulter le Mode d'Emploi
 Gebrauchsanweisung
 Consultare le istruzioni per l'uso
 Consultar el manual de instrucciones
 Consultar o Folheto de Instruções
 Se brugsanvisningen

In Vitro Diagnostic Medical Device
 Produit Diagnostique Médical In Vitro
 In-vitro-Diagnostika
 Diagnostico in vitro Dispositivo medico
 Producto Sanitario para Diagnóstico In Vitro
 Produto de Diagnóstico Médico In vitro
 Medicinsk udstyr til in vitro-diagnostik

Temperature Limitation
 Températures Limites
 Temperaturbegrenzung
 Temperatura limite
 Limite de temperatura
 Temperatura Limitada
 Opbevarings ved

Biological Risks
 Risques Biologiques
 Biogefährdung
 Rischi biologici
 Riesgo biológico
 Riscos Biológicos
 Biologisk risiko

Manufacturer
 Fabricant
 Hersteller
 Fabricante
 Fabricante
 Fabricante
 Fabrikant

Cellabs Pty Ltd
 Unit 7, 27 Dale Street (PO Box 421)
 Brookvale, NSW 2100 Australia
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426
 Web: http://www.cellabs.com.au
 Email: sales@cellabs.com.au

Authorised Representative in the EC
 Représentant Autorisé dans la EC
 Bevollmächtigter in der Europäischen Gemeinschaft
 Rappresentante Autorizzato nella EC
 Representante autorizado en la EC
 Representante Autorizado na CE
 Repräsentant

MDCI Ltd
 Arundel House, 1 Liverpool Gardens
 Worthing, West Sussex BN11 1SL
 United Kingdom

Insert Version
 Manuel Version
 Fassung der Packungsbeilage
 Versione Inserto
 Versión del manual de instrucciones
 Versão do folheto de instruções
 Indlægsæddel, version

LR2.10 -PDF
 28th November 2006

INTENDED USE AND PRINCIPLE OF THE TEST

The Crypto/Giardia Cel IF Test is an *in vitro* direct immunofluorescence test for the simultaneous detection of *Cryptosporidium* oocysts and *Giardia* cysts in faecal and environmental specimens. The fluorescein-labelled mouse monoclonal antibody reagent binds specifically to *Cryptosporidium* oocysts and/or *Giardia* cysts in the specimen. Organisms display bright green fluorescence with typical morphology.

CONTENTS OF THE KIT

	RR2	KR2 Demo	KR2 Standard	Bulk
	Crypto/Giardia Cel Reagent	0.25mL	1.25mL	5mL
	Positive Control Slide (Single Use only)	1	1	-
	Mounting Fluid	2.5mL	2.5mL	-
	<i>Tests</i>	10	50	200

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; acetone for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x200-x1000 magnification.

PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix components from different kits. The Crypto/Giardia Cel Reagent has been optimised for use with Cellabs Positive Control Slide and Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. The detection of *Cryptosporidium* and *Giardia* from water is dependent on the efficient filtration and recovery from large volumes of water (500-1000L). A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

SPECIMEN PREPARATION

Stool specimens, either fresh or preserved in 10% formalin or SAF may be used (Specimens treated with other fixatives have not been tested). A procedure for processing fresh stool specimens is described below. If stool specimens are not processed immediately they may be stored at 2-8°C for 24-48 hours. Stools may be concentrated prior to processing using concentration kits. They may also be concentrated by the formalin-ether sedimentation method.

Prepare an approximate 1/10-1/50 dilution of stool. Dilute 50µL or 5mm diameter of stool in PBS with 0.1% sodium azide and mix thoroughly to disperse the specimen using an applicator stick and allow large stool particles to sediment.

PREPARATION OF SLIDES

- Place 20µL of the faecal or concentrated specimen from water onto a microscope slide.
- Allow the specimen to completely air dry.
- Fix the slides for 5 minutes in acetone and allow to air dry.

INSTRUCTIONS FOR USE

- Add 25µL of RR2 to the fixed specimen and CG, covering well area.
- Incubate the slides at 37°C in a humid chamber for 30 minutes. Do not allow the slides to dry as this may cause non-specific binding.
- Rinse gently in a bath of PBS for one minute.
- Drain slide and remove excess moisture around well with tissue.
- Add a drop of RMG to the slide well. Place a coverslip on top of the drop and remove air bubbles.
- Scan the entire specimen using a fluorescence microscope initially at X200 magnification, then at X400, and X1000 for confirmation. Read immediately or store at 2-8°C in the dark for up to 24 hours.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

Cryptosporidium oocysts are 2-6µm in size, round or oval in shape with bright green fluorescence. A fold or suture can be seen on the surface. This suture is less apparent in environmental samples. *Giardia* cysts are 8-12µm in size, elliptical in shape with bright green fluorescence. The test may be considered positive if one or more oocysts or cysts are present. The organisms, if present together in the same specimen, may be distinguished by the difference in their size and shape. Both organisms are seen against a background of reddish brown counterstained material. A comparison should be made with the positive control slide for the appearance of the organisms. Negative specimens should have signs of faecal material containing debris such as undigested food particles, cells (counterstained red-brown), and possibly other organisms such as bacteria and protozoa, and should not show any organisms of characteristic green fluorescence and morphology. Non-specific fluorescence may occur in some specimens but this may be distinguished morphologically from the *Cryptosporidium* oocysts and *Giardia* cysts. Algae may be present in environmental specimens and these fluoresce red-orange.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE CRYPTO/GIARDIA CEL

Refer to summary table at end of insert. All data on the Crypto/Giardia Cel can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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