

**FIGURE 1: POSITIVE CONTROL SLIDE & MOUNTING FLUID DIAGRAM FOR USE**

Please refer to individual Cellabs IFA Kit Insert when performing test

**For Chlamydia Cel Pn IFA insert the following steps**

**EXPLANATION OF SYMBOLS**

Consult Instructions for Use  
 Consulter le Mode d'Emploi  
 Gebrauchsanweisung  
 Consultare le istruzioni per l'uso  
 Consultar el manual de instrucciones  
 Consultar o Folheto de Instruções  
 Se brugsanvisningen  
 Συμβουλευτείτε τις Οδηγίες Χρήσης  
 Tanulmányozza a használati utasítást  
 Zapoznać się z instrukcją użycia

In Vitro Diagnostic Medical Device  
 Produit Diagnostique Médical In Vitro  
 In-vitro-Diagnostika  
 Diagnostico in vitro Dispositivo medico  
 Producto Sanitario para Diagnóstico In Vitro  
 Produto de Diagnóstico Médico In vitro  
 Medicinsk udstyr til *in vitro*-diagnostik  
 In Vitro Διαγνωστική Ιατρική Συσκευή  
 In Vitro Orvosi Diagnosztikai Eszköz  
 Wyrób medyczny do diagnostyki in vitro

Temperature Limitation  
 Températures Limites  
 Temperaturbegrenzung  
 Temperatura limite  
 Limite de temperatura  
 Temperatura Limitada  
 Opbevaras ved  
 Περιορισμός Θερμοκρασίας  
 Tárolási hőmérséklet  
 Ograniczenia temperatury

Biological Risks  
 Risques Biologiques  
 Biogefährdung  
 Rischi biologici  
 Riesgo biológico  
 Riscos Biológicos  
 Biologisk risiko  
 Βιολογικοί Κίνδυνοι  
 Biológiai veszély  
 Ryzyko biologiczne

Authorised Representative in the EC  
 Représentant Autorisé dans la EC  
 Bevollmächtigter in der Europäischen Gemeinschaft  
 Rappresentante Autorizzato nella EC  
 Representante autorizado en la EC  
 Representante Autorizado na CE  
 Représentant  
 Εξουσιοδοτημένος Αντιπρόσωπος στην EC  
 Felhatalmazott képviselő az EU-ban  
 Autoryzowany przedstawiciel w UE


MDCI Ltd  
 Arundel House, 1 Liverpool Gardens  
 Worthing, West Sussex BN11 1SL  
 United Kingdom

Manufacturer  
 Fabricant  
 Hersteller  
 Fabbricante  
 Fabricante  
 Fabricante  
 Fabrikant  
 Κατασκευαστής  
 Gyártó  
 Producent

Cellabs Pty Ltd  
 Unit 7, 27 Dale Street (PO Box 421)  
 Brookvale, NSW 2100 Australia  
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426  
 Web: <http://www.cellabs.com.au>  
 Email: [sales@cellabs.com.au](mailto:sales@cellabs.com.au)

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LX.04 - PDF  
 28th November 2006

**POSITIVE CONTROL SLIDE & MOUNTING FLUID**

The following Instructions for Use provides information on Cellabs' Positive Control Slides and Mounting Fluids included in Cellabs IFA tests and available individually

**INTENDED USE AND PRINCIPLE OF THE DEVICES**

The Positive Control Slide is to be used as a positive control when performing an *in vitro* direct immunofluorescence test (using Cellabs' IF Reagent) for the detection and diagnosis of tropical and infectious organisms in patient specimens and environmental samples. Organisms display bright green fluorescence with typical morphology. The Mounting Fluid is used to mount coverslip to specimen slide when viewing under a fluorescent microscope. Refer to Table 1 for intended use of Positive Control Slide and Mounting Fluid.

**PACKAGING**

Positive Control Slide is individually packaged in foil bag with one desiccant sachet, and is ready for use. Mounting Fluid is packaged in plastic vials containing 2.5mL or 10mL and is ready for use. Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and do not change once opened.

**MATERIALS THAT MAY BE REQUIRED BUT ARE NOT PROVIDED**

Microscope slides with 6-8mm diameter wells; acetone or methanol for specimen fixation (refer to Kit insert or Reagent insert for which fixative to use); Cellabs' IF Reagent; precision pipette for delivering 25µL; humid chamber; wash bath; (PBS) for washing step; cover slips; non-fluorescing immersion oil; fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x10-x1000 magnification.

**PRECAUTIONS**

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Patient specimens should be handled as though potentially infectious. The positive control slide has been fixed and therefore is non-viable, however should always be handled as though potentially infectious. A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

**SPECIMEN COLLECTION AND PREPARATION OF SPECIMEN SLIDES**

Refer to Kit Package Insert for detailed instructions on using Positive Control Slide and Mounting Fluid with Cellabs IF Reagent.

**INSTRUCTIONS FOR USE**

Refer to Kit Package Insert for detailed instructions on using Positive Control Slide and Mounting Fluid with Cellabs IF Reagent.

**READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS**

Refer to Kit Package Insert for reading and interpreting results of using Positive Control Slide and Mounting Fluid with Cellabs IF Reagent.

**WASTE DISPOSAL**

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

**SENSITIVITY, SPECIFICITY, & OTHER DATA ON POSITIVE CONTROL SLIDE AND MOUNTING FLUID**

Not applicable.

**INDEMNITY NOTICE**

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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