**International Symbols**  
**Symboles internationaux**  
**Internationale symbolen**  
**Simboli internazionali**  
**Simbolos internacionales**  
**Simbolos Internacionais**  
**Internationale symboler**  
**Διεθνή Σύμβολα**  
**Nemzetközi szimbólumok**

Consult Instructions for Use  
*Consulter le mode d'emploi*  
*Raadpleeg de gebruiksaanwijzing*  
*Consultare le istruzioni per l'uso*  
*Consultar instrucciones de uso*  
*Consulte as Instruções de Utilização*  
*Se brisingenugsanv*  
*Συμβουλευτείτε τις Οδηγίες Χρήσης*  
*Tekintse meg a használati utasítást*

Do Not Re-use  
*Ne pas réutiliser*  
*Niet hergebruiken*  
*Non riutilizzare*  
*No reutilizar*  
*Não reutilizar*  
*Må ikke genbruges*  
*Μην το ξαναχρησιμοποιήσετε*  
*Ne használja újra*

In Vitro Diagnostic Medical Device  
*Dispositif médical de diagnostic in vitro*  
*Medisch diagnostisch hulpmiddel voor in vitro gebruik*  
*Dispositivo medico diagnostico in vitro*  
*Dispositivo médico para diagnóstico in vitro*  
*Inserir versão Dispositivo Médico para Diagnóstico In Vitro*  
*In vitro-diagnostisk medicinsk udstyr*  
*Διαγνωστική Ιατρική Συσκευή In Vitro*  
*In vitro diagnosztikai orvostechnikai eszköz*

Use By/Expiration Date  
*Date de péremption*  
*Vervaldatum*  
*Data di scadenza*  
*Fecha de caducidad*  
*Data de Validade*  
*Udløbsdato*  
*Ημερομηνία Λήξης*  
*Lejárati dátum*

Temperature Limitation  
*Limitation de temperature*  
*Temperatuurlimiet*  
*Limite di temperature*  
*Limite de temperatura*  
*Limitação de Temperatura*  
*Temperaturbegrænsning*  
*Περιορισμός Θερμοκρασίας*  
*Hőmérsékletkorlát*

Cellabs Pty Ltd  
 Unit 7, 27 Dale Street  
 Brookvale, NSW 2100 Australia  
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426  
 Web: <http://www.cellabs.com.au>  
 Email: [sales@cellabs.com.au](mailto:sales@cellabs.com.au)


Batch  
 Lot  
 Batch  
 Lotto  
 Lote  
 Lote  
 Batch  
 Παρτίδα  
 Tétel

WMDE B.V  
 Bergerweg 18  
 6085 AT Horn  
 The Netherlands

Control Positive  
*Contrôle positif*  
*Positieve controle*  
*Controllo positivo*  
*Control positivo*  
*Controlle Positivo*  
*Positiv control*  
*Θετικός Έλεγχος*  
*Pozitív kontroll*

IFU Version  
 IFU la version  
 IFU Versie  
 IFU la versione  
 Versión insertada  
 IFU versão  
 IFU version  
 Εισαγάγετε την έκδοση  
 IFU verzió

LC1.17  
 21 January 2026

**CHLAMYDIA CEL**

INSTRUCTIONS FOR USE

	English	2
	Français	3
	Deutsch	4
	Italiano	5
	Español	6
	Português	7
	Dansk	8
	Ελληνικά	9
	Magyar	10

Sensitivity, Specificity & Other Data	11
Explanation of Symbols	12

*Product Code: KC1*

# CHLAMYDIA CEL

## INTENDED USE AND PRINCIPLE OF THE TEST

The Chlamydia Cel IFA is a direct immunofluorescence assay (IFA) kit. It is a rapid, qualitative, in vitro test for detecting *Chlamydia trachomatis* organisms in clinical specimens. The fluorescein-labelled mouse monoclonal antibody reagent binds specifically to the major outer membrane protein (MOMP) of *C. trachomatis*. It can be used for specimens prepared as direct patient smears or tissue culture confirmation samples fixed on diagnostic slides. This test is not automated and is intended for use by laboratory professionals or personnel trained and competent in immunodiagnostic methods.

## CONTENTS OF THE KIT

		KC1	Bulk
	Chlamydia Cel Reagent	1.25mL	5mL
	Positive Control Slide (Single Use only)	1	1
	Mounting Fluid	2.5mL	10.0 mL
	Tests	50	200

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

## MATERIALS REQUIRED BUT NOT PROVIDED

Microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; methanol for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x600-x1000 magnification.

## PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix components from different kits. The Chlamydia Cel Reagent has been optimised for use with Cellabs Positive Control Slide and Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information. For professional use only. To be used with confirmatory testing guidelines or standard quality assurance systems for clinical samples.

## SPECIMEN COLLECTION

**Cervical Specimens:** Specimens from the female cervix should contain as many columnar epithelial cells as possible as *C. trachomatis* are intracellular organisms that infect these cells. Clean cervix with a sterile cotton swab before sampling. Insert swab in the cervical canal and rotate swab at the squamo-columnar junction. Withdraw swab without touching any vaginal surfaces.

**Urethral Specimens:** Specimens from the male urethra should also contain intact epithelial cells to ensure an accurate diagnosis. Patients preferably should not have urinated one hour before collection. Insert swab 2-4cm into the urethra using dacron or cotton wool-tipped thin swabs. Rotate the swab and withdraw.

**Conjunctival Specimens:** Using a sterile cotton-tipped swab, vigorously scrape the lower lid conjunctivae of affected eyes to obtain a good yield of epithelial cells. If both eyes are affected, first swab the least affected eye then swab the most affected eye.



**Culture Confirmation:** Immerse the swab in chlamydia transport medium and use for inoculating coverslip or monolayer cultures using standard techniques.

## PREPARATION OF SLIDES

A. Direct Specimens: Roll the swab back and forth on the diagnostic glass slide to cover the area of the well. Allow the specimen to air dry, then fix in methanol for 5 minutes. Drain and air dry. If the specimen is not tested immediately, store at 2-8°C overnight or freeze at -20°C for up to 2 months.

B. For culture specimens: Fix coverslip or monolayer culture sample on slide with methanol for 5 minutes. Drain and dry ready for staining.

## INSTRUCTIONS FOR USE

- Add 25µL of  to the fixed specimen and positive control slide , covering well area.
- Incubate the slides at 37°C in a humid chamber for 30 minutes. Do not allow the slides to dry as this may cause non-specific binding.
- Rinse gently in a bath of PBS for one minute.
- Drain slide and remove excess moisture around well with tissue.
- Add a drop of  to the slide well. Place a coverslip on top of the drop and remove air bubbles.
- Scan the entire specimen using a fluorescence microscope under oil immersion at x600-x1000 magnification. Read immediately or store at 2-8°C in the dark for up to 24 hours.

## READING AND INTERPRETATION OF RESULTS

Extracellular elementary bodies (EB) appear as bright apple-green fluorescent pin-point, smooth-edged disc shaped bodies (approximately 300nm in diameter) and can be seen against a background of reddish-brown counterstained cells. Reticulate bodies may also be observed and are 2-3 times larger than the EB's and they fluoresce evenly or possess dark centres with a halo of fluorescence. Intact intracellular chlamydial inclusions are rarely seen. The control slide should be used for comparison with the appearance and size of EB's found in the specimen. A positive diagnosis is made when fixed stained specimens show at least 10 chlamydial bodies. A negative diagnosis is reported when fixed stain smears are free of chlamydial organisms but at least 10 columnar cells (intact or ruptured) are present. Irregularly shaped fluorescent material that differs in size from chlamydial bodies described above or fluoresces white, red or yellow should be considered non-specific staining.

## WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

## SENSITIVITY, SPECIFICITY, & OTHER DATA ON CHLAMYDIA CEL

Refer to summary table at end of insert. All data on the Chlamydia Cel can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

## INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

## TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE CHLAMYDIA CEL

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