CHLAMYDIA CEL LPS

English
Product Code: KC2

FIGURE 1: CHLAMYDIA CEL LPS DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE CHLAMYDIA CEL LPS

TABELEAU 1: SENSIBILITÉ, SPÉCIFICITÉ ET AUTRES DONNÉES DU TEST CHLAMYDIA CEL LPS
 TABELLE 1: SENSIBILITÄT, SPEZIFITÄT UND ANDERE DATEN ZUM CHLAMYDIA CEL LPS
 TABELLA 1: SENSIBILITÀ, SPECIFICITÀ ED ALTRI DATI SULLA CHLAMYDIA CEL LPS
 TABLA 1: SENSIBILIDAD, ESPECIFICIDAD Y OTROS DATOS DEL CHLAMYDIA CEL LPS
 TABELA 1: SENSIBILIDADE, ESPECIFICIDADE, & OUTROS DADOS DO CHLAMYDIA CEL LPS

Trial Essai Versuch Prova Prueba Teste	Sensitivity Sensibilité Sensitivität Sensibilità Sensibilidad Sensibilidade	Specificity Spécificité Spezifität Specificità Especificidad Especificidade	Repeatability Répétabilité Wiederholpräzision Ripetibilità Repetibilidad Repetição	Reproducibility Reproductibilité Reproduzierbarkeit Riproducibilità Reproducibilidad Reprodutibilidade
A	89%	94%	-	-
B	88%	100%	-	-
C	90%	84%	-	-
D	-	-	100% Correlation	100% Correlation
<p>Not cross reactive with / Pas de Réaction Croisée avec / Keine Kreuzreaktionen mit / Non mostra reazione crociata con / No muestra reacción cruzada con / Sem reacção cruzada com:</p> <p><i>Staphylococcus epidermidis, Staphylococcus aureus, Beta haemolytic Streptococcus group B, Beta haemolytic Streptococcus group A, Lactobacillus acidophilus, Peptostreptococcus anaerobius, Proteus mirabilis, Streptococcus faecalis, Streptococcus pneumoniae, Haemophilus influenzae, Haemophilus parainfluenzae, Escherichia coli, Klebsiella pneumoniae, Proteus mirabilis, Pseudomonas aeruginosa, Bacteriodes fragilis, Viridans Streptococcus, Candida albicans, Bacteriodes melaninogenicus, Actinomyces israelii, Gardnerella vaginalis, Neisseria gonorrhoeae, Trichomonas vaginalis:</i></p>				

EXPLANATION OF SYMBOLS

- Consult Instructions for Use
- In Vitro Diagnostic Medical Device
- Temperature Limitation
- Batch
- Control Positive
- Use By/Expiration Date
- Do Not Re-use

Cellabs Pty Ltd
 Unit 7, 27 Dale Street
 Brookvale, NSW 2100 Australia
 Tel: +61 2 9905 0133 Fax: +61 2 9905 6426
 Web: <http://www.cellabs.com.au>
 Email: sales@cellabs.com.au

Insert
 Version

en fr de it es pt
 LC2.10.1
 23 May 2017

INTENDED USE AND PRINCIPLE OF THE TEST

The Chlamydia Cel LPS IF Test is a rapid *in vitro* direct immunofluorescence test for the detection of *Chlamydia* organisms in a variety of specimens. The fluorescein-labelled mouse monoclonal antibody reagent binds specifically to *Chlamydia* lipopolysaccharide (LPS) in the specimen. Chlamydial organisms display bright green fluorescence with typical morphology.

CONTENTS OF THE KIT

		KC2 Kit	Bulk
	Chlamydia Cel LPS Reagent	1.25mL	5mL
	Positive Control Slide (Single Use only)	1	-
	Mounting Fluid	2.5mL	-
	Tests	50	200

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; acetone or methanol for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x600-x1000 magnification.

PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix components from different kits. Chlamydia Cel LPS Reagent has been optimised for use with Cellabs Positive Control Slide and Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

SPECIMEN COLLECTION

Samples should be collected by swabbing using cotton tipped swabs and should contain as many epithelial cells as possible. Purulent discharge should be cleaned from the site before sampling.

PREPARATION OF SLIDES

Swabs should be applied to a microscope slide by rolling back and forth over a small area (approximately 10mm in diameter). Allow specimens to air dry and fix in acetone or methanol for 5 minutes. Drain and air dry. If the specimen is not tested immediately, store at 2-8°C overnight or freeze at -20°C for up to 2 months.

INSTRUCTIONS FOR USE

1. Add 25µL of to the fixed specimen and , covering well area.
2. Incubate the slides at 37°C in a humid chamber for 30 minutes. Do not allow the slides to dry as this may cause non-specific binding.
3. Rinse gently in a bath of PBS for one minute.
4. Drain slide and remove excess moisture around well with tissue.
5. Add a drop of to the slide well. Place a coverslip on top of the drop and remove air bubbles.
6. Scan the entire specimen using a fluorescence microscope under oil immersion at x600-x1000 magnification. Read immediately or store at 2-8°C in the dark for up to 24 hours.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

The most common chlamydial forms seen in specimens are extra-cellular elementary bodies. They appear as bright apple-green fluorescent pinpoint, smooth-edged disc shaped bodies (approximately 300nm in diameter) and can be seen against a background of reddish-brown counterstained cells. Reticulate bodies may also be observed. These are 2-3 times larger than the EBs and they either fluoresce evenly or possess dark centres with a halo of fluorescence. Intact intracellular Chlamydia inclusions will be seen in eye specimens and impression smears. Any material which can be distinguished from chlamydial forms or fluoresce other than apple-green should be disregarded. The control slide should be used for comparison with the appearance and size of EB's found in the specimen.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON CHLAMYDIA CEL LPS

Refer to summary table at end of insert. All data on Chlamydia Cel LPS can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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