# ENTAMOEBA CELISA PATH

## INTENDED USE AND PRINCIPLE OF THE TEST

The Entamoeba CELISA Path test is an enzyme immunoassay for the rapid detection of the adhesin of *E. histolytica* in human faecal specimens. The Entamoeba CELISA Path test uses antibodies to the adhesin of the organism. The microtitre wells in the kit contain immobilised polyclonal antibody that binds adhesin. The conjugate is a monoclonal antibody specific for the adhesin of *E. histolytica*. In the assay an aliquot of a faecal specimen is emulsified in the diluent and the diluted specimen is transferred to a microtitre well. If adhesin is present in the specimen, it binds to the conjugate and immobilised polyclonal antibody during the incubation phase. Any unbound material is removed during the washing steps. Following the addition of substrate, a colour is detected due to the enzyme-antibody-antigen complexes that form in the presence of adhesin.

## CONTENTS OF THE KIT

<b>EPMW</b>	Celisa Plate – 1x96 wells - (single use only)		1 plate
<b>CONTROL</b>	Positive Control	(Black Cap)	3.5mL
<b>EPDL</b>	Diluent	(White Cap)	40mL
<b>EPPO</b>	Enzyme Conjugate	(Red Cap)	7.0mL
<b>EPWB</b>	Wash Buffer Concentrate (20x)	(White Cap)	50mL
<b>EPSA</b>	Substrate	(Blue Cap)	14.0mL
<b>EPSS</b>	Stopping Solution	(Yellow Cap)	7.0mL
	Plastic Adhesive Sheets		1

Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

## MATERIALS REQUIRED BUT NOT PROVIDED

Distilled water for diluting wash reagent; 1L bottle for diluting wash reagent; Squirt bottle for wash reagent; Test tubes; Applicator sticks; Micropipettes with disposable tips; Vortex mixer; Timer; Discard container and absorbent paper; Spectrophotometer capable of reading absorbance of EIA plates at 450nm or 450/620nm.

## PRECAUTIONS

For in vitro diagnostic use only. Reagents should not be used after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Caps and tips are colour coded, do not mix. When handling assay wells; avoid scratching the bottom of the wells as this may result in elevated absorbance readings. Thimerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eyes and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. The positive control has been shown to be non-infectious in tissue culture; however it should be handled as though potentially infectious. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Material Safety Data Sheet.

## INSTRUCTIONS FOR USE

### Preparation of Wash Buffer

If crystals are present in the concentrate, warm to dissolve. Prepare working strength **WASH BUFFER** by adding 50mL of **EPWB** to 950mL of distilled water. Label the bottle and store at 2-8°C.

### Specimen Collection and Preparation

Standard stool specimen collection and handling procedures for culture are appropriate. **Do not use specimens that have been collected or stored in 10% formalin, SAF, or PVA fixatives.** Stool specimens that are in transport media such as Cary Blair or C&S may be used; however, these samples arrive at a 1:5 dilution and must be assayed straight from the container and not diluted in **EPDL**. Specimens should be transported as soon as possible and stored at 2 to 8°C. If possible, test stool specimens which are less than 24 hours old. If the test cannot be performed within 48 hours, store specimens at -20°C or lower. Freezing and thawing of the specimen, especially multiple times, may result in loss of activity due to degradation of the adhesin.

Make sure that specimens are thoroughly mixed prior to performing the assay. This includes complete mixing of the specimen prior to transfer to diluent, as well as complete mixing of the diluted specimen prior to transfer to the microwell. The diluent has been formulated to stabilise the adhesin in faecal specimens and minimise degradation. When test is ready to be performed, dilute samples as follows:

1. Set up one vial for each specimen to be tested. Add 0.4mL **EPDL** to each vial. Label the vial directly on the side.
2. Thoroughly mix (vortex) the faecal specimen to ensure adequate sampling.
3. **Formed stools:** use a swab to transfer the faecal specimen to the vial. Coat the swab completely before transferring the specimen. Mix the swab in the **EPDL** to remove as much sample as possible and squeeze the swab against the side of the vial to express any residual liquid. The procedure results in the transfer of approximately 0.15 to 0.20g of specimen. If the swab is not used for formed stools, transfer approximately 0.15 to 0.20g of specimen to the **EPDL**.

**Liquid stools:** use a plastic pipette to transfer 0.4ml specimen to vial. Make sure the liquid specimens are evenly suspended before transferring.

4. Vortex the vials for 10 seconds and store at 2 to 8°C until the ELISA is performed. Diluted specimens (**SAMPLE**) should be used within 5 hours of preparation. Vortex again before transferring **SAMPLE** to microtitre well. This ensures thorough mixing of the specimen.

### Assay Procedure

1. Bring all reagents to room temperature (18-25 °C) before use.
2. Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
3. Remove required number of **EPMW** strips. Reseal the foil bag containing unused microwell strips immediately with tape.
4. Add 1 drop of **EPPO** (red cap) to a positive control well, negative control well and patient sample well. Hold the conjugate vial vertically when adding the drops. Identification marks may be written directly on side of well.
5. Two (2) control wells must be used each time the test is performed. Add one (1) drop of the positive control **CONTROL** to the positive control well and 0.1mL of diluent **EPDL** to the negative control well. Add 0.2mL of **SAMPLE** to the test well(s). Change pipette tips between all samples and controls. Cover the wells and incubate at room temperature for two (2) hours.
6. Shake out the contents of the assay wells. Wash the wells by using the **WASH BUFFER** in a squirt bottle with fine-tipped nozzle, directing the wash buffer to the bottom of the well with force. Fill the wells, and then shake the wash buffer out. Bang the inverted plate on a dry paper towel and repeat the washing step another three (3) times. If any particulate matter is in the wells continue to wash until it is removed.
7. After washing, bang dry (until no liquid comes out) and dispose of paper towels and specimen containers.
8. Dispense two (2) drops of substrate **EPSA** (blue cap) to each well. Gently tap the wells to mix the substrate. Incubate for ten (10) minutes at room temperature in the dark. Tap the wells two (2) or three (3) times during incubation to mix the contents. Observe the blue colour formation.
9. Dispense one (1) drop of stopping solution **EPSS** (yellow cap) to each well. The addition of the stop solution changes the colour of the wells from blue to yellow.
10. Quantitate the results by measuring the absorbance in a spectrophotometer at a wavelength of 450nm (reference 620nm on dual wavelength reader). The instrument should be blanked against air. Wipe the underside of each well before measuring the optical density. Plates should be read within 10 minutes of stopping the reaction.

If a spectrophotometer is unavailable, the test may be read visually in good light against a white background (such as a sheet of paper)

## READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

### Visually

**Negative:** Any sample that is not obviously more yellow than the negative control well.

**Positive:** Any sample that is obviously more yellow than the negative control well.

**Note:** The negative control, as well as some test wells, may show some slight yellow colour. A sample well must be obviously darker than the negative control well to be called a positive result.

### Photometrically

	450nm	450/620nm
Negative Sample & Negative Control	OD <sub>450</sub> < 0.150	OD <sub>450/620</sub> < 0.090
Positive Sample	OD <sub>450</sub> ≥ 0.150	OD <sub>450/620</sub> ≥ 0.090
Positive Control (once negative control value has been subtracted)	OD <sub>450</sub> ≥ 0.50	

Values for positive and negative control must be within range for the test to be valid.

A positive test result indicates that *E.histolytica* adhesin is present in the faecal specimen. A negative result indicates that *E.histolytica* adhesin is not present in the faecal specimen.

### WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

### SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE ENTAMOEBIA CELISA PATH

Refer to summary table at end of insert. All data on the Entamoeba CELISA Path can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

### INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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Cellabs Pty Ltd