FILARIASIS ANTIGEN Og4C3 ELISA

ELISA Kit for the detection and quantification of *Wuchereria bancrofti* Antigen

English

INTRODUCTION

Lymphatic filariasis in man is caused by infection with the filarial parasites *W. bancrofti* and *Brugia* spp. These parasites inhabit the lymphatics and cause disease by obstructing and inducing secondary inflammatory changes in lymph vessels. The disease is transmitted by mosquitoes, which ingest microfilariae during feeding on an infected host and transmit the infective larvae to other individuals at a subsequent feeding. Worldwide, more than 90 million people are affected by lymphatic filariasis. Most live in the humid tropics in areas such as Africa (south of the Sahara), Egypt, the Indian subcontinent, South-East Asia, China, Madagascar, Papua New Guinea, the Pacific Islands, the Philippines and Central and South America (World Health Organisation, 1984).

Acute symptoms of lymphatic filariasis primarily involve lymphadenitis and lymphangitis. Recurrent fever and pain in affected lymph nodes are the usual sequelae. In some patients, symptoms may be less specific, with only fever and malaise present. Some infected individuals are asymptomatic. Lymphoedema becomes apparent after repeated episodes of lymphadenitis, and swelling of the limbs or scrotum may occur. Lymphoedema and elephantiasis may affect the leg, arm and scrotum and occasionally the vulva and breasts, but they are usually restricted to the leg below the knee. With brugian filariasis, more severe inflammatory changes are noted in the lymphatics, whereas bancroftian filariasis has a more extensive swelling of the entire limb.

Filariasis in previously unexposed migrants to an endemic area has a similar clinical course, but can be manifested earlier (6 - 8 weeks) than the normal 7 - 8 month clinical incubation period. Microfilaraemia is an uncommon finding in these migrant individuals, reflecting an intense immunological reaction to the parasite (Partono et al., 1977).

INTENDED USE

The TropBio Filariasis Antigen Og4C3 ELISA was developed for the detection of *Wuchereria bancrofti* (Bancroftian filariasis) circulating filarial antigen (CFA) using whole blood, plasma, serum or dried blood samples (DBS).

This assay can be used as a diagnostic tool for clinical cases, screening or for epidemiological surveys to determine the prevalence of disease in study sites based on the detection of *W. bancrofti* CFA in populations. The assay is increasingly used as a tool to assess and monitor the prevalence of disease after rounds of mass drug administration (MDA). In low prevalence settings, the interpretation of ELISA result can be modified to focus on the low antigen readings and define an appropriate cut off value as an alternative to using the "titre group method" by Moore and Copeman (1990).

CONTENTS OF THE KIT

	Code	KF1 (1 plate)	KF2 (5 plates)
Microwell plates (Og4C3)	FGMW	1	5
Sample Diluent (1x)	FGSD	30 mL	120 mL
Antibody and Conjugate Diluent (1x)	FGCD	30 mL	100 mL
Standard Antigens (1-7)	FGPC	-	0.8 mL
Standard Antigen No. 7 (Positive Control)	FGPC	0.6 mL	-
Standard Antigen No. 1 (Negative Control)	FGNC	0.6 mL	-
Rabbit anti-onchocerca antibody (120x)	FGAB	0.07 mL	0.35 mL
Anti-rabbit HRP-conjugate (120x)	FGPO	0.07 mL	0.35 mL
TMB Substrate Concentrate (20x)	FGSC	0.5 mL	1.8 mL
Substrate Buffer (1x)	FGSB	7.0 mL	30.0 mL
Stopping Solution (1x)	FGSS	7.0 mL	30.0 mL
Wash Buffer (20x)	FGPT	100 mL	250 mL

Store kit at 2-8°C. Store the positive control frozen using the instructions below. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Distilled water, micropipettes and tips, clean glassware or plastic containers for solutions, humid chamber, ELISA washer, ELISA plate reader to read absorbances at a single wavelength of 450nm or at dual wavelength of 450nm and 620nm or automated ELISA systems.

PRECAUTIONS

For in-vitro diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different batch numbers. Thiomerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eye and mucous membranes. Ensure all dilution containers are clean before use. Avoid exposure of the substrate to light. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Safety Data Sheet. Reagents should not be used after the expiry date shown on the label.

INSTRUCTIONS FOR USE

METHOD FOR SERUM OR PLASMA SAMPLES

All steps carried out at room temperature. Ensure that all reagents and the microwell plates are at room temperature before use. **Quantities indicated below refer to the use of ONE PLATE. Calculate the required reagent depending on the wells being used.** To prepare wash buffer, dilute 25mL of Wash Buffer 20x concentrate to 475mL distilled water. Wash Buffer 1x can be stored at 2-8°C for up to one week.

The boiling treatment is optional. See Notes on the Assay.

Sample Preparation

To prepare samples (serum or plasma), add 100µL of test sample to 300µL of sample diluent in a suitable tube for boiling eg. Eppendorf microcentrifuge tubes (piece the lid with a fine needle to allow air to escape during boiling). Do not boil or dilute standards 1-7.

Place the tubes into a 100°C boiling water bath for five minutes. After boiling, centrifuge the samples at 2,000g for 15 minutes (racked tubes) or 10,000g for five minutes (Eppendorf tubes). The clear supernatant fluid contains the heat stable antigen.

1. Add 50µL aliquots of boiled sample supernatant fluid to a test well. Up to 80 samples can be tested per plate.
Add 50µL of Standard Antigens (1-7) and conjugate control (use sample diluent) in duplicates using two microwell strips. Refer to Figure 1 for the plate layout diagram.
2. Incubate by placing the plate in a humid container for 60 minutes at 37°C. (Note: using an overnight incubation will result in high OD readings that will not produce a typical standard curve).
3. Wash the plate three times with wash buffer, invert and tap gently to remove residual buffer.

4. Dilute 50uL of the anti-Onchocerca antibody (Yellow cap) to 6 mL of antibody diluent (Blue solution). Mix thoroughly. Add 50 uL of diluted rabbit-anti Onchocerca antibody to all wells. Place the plate in a humidity chamber and incubate for 45 minutes at 37°C .
6. Wash the plate three times as before.
7. Dilute 50 uL of Anti-rabbit HRPO conjugate (Purple cap) to 6 mL of antibody diluent. Mix thoroughly. Add 50 uL of diluted conjugate to all wells. Place in a humidity chamber and incubate for 45 minutes at 37°C.
8. Wash the plate three times as before.
9. Prepare the substrate by diluting 300uL of TBM substrate concentrate to 5700uL of Substrate Buffer. Mix thoroughly. Add 50uL to each well and incubate for 15 minutes in the dark at room temperature.
10. Stop the reaction by adding 50uL of Stopping Solution to each well (blue colour will turn yellow).
11. Read the plate using a plate reader at 450 nm, or dual wavelengths of 450/620 nm.

METHOD FOR FILTER PAPER SAMPLES

Note: This kit does not contain 1% hydrogen peroxide (H2O2). You are required to prepare the solution in Step 5 of the method. Filter paper discs can be ordered separately. Note that using different types of filter paper other than the Tropbio filter paper discs may give different results due to various paper thickness, size, make and pre-treatments.

1. **Preparation of filter paper samples:** Cut three protrusions from the filter paper disk and add them to a suitable tubes eg. 2 mL Eppendorf microtubes. Protrusions can be cut in half to allow the disks to reach the bottom of the tubes if not using 2mL tubes.
Add 200 uL of sample diluent to each tube and leave the blood samples to elute overnight at 2-8°C. The next morning, place the tubes into a 100°C water bath for five minutes
After heating in the water bath, centrifuge the samples at 2,000 g for 15 minutes. The supernatant contains the heat stable antigen.
2. Add 50uL aliquots of boiled sample supernatant fluid to a test well. Up to 80 samples can be tested per plate.
Add 50uL of Standard Antigens 1-7 (do not boil) and conjugate control (use sample diluent) in duplicates using two strips of a plate. Refer to the plate layout diagram on Page 4.
3. Incubate by placing the plate in a humid container for 60 minutes at 37°C . (*Note: using an overnight incubation will result in high OD readings that will not produce a typical standard curve*).
4. Wash the plate three times with wash buffer, invert and tap gently to remove residual buffer.
5. Prepare a 1% hydrogen peroxide solution for ONE plate by adding 400uL of hydrogen peroxide (~30%) to 12mL of 1x wash buffer.
6. Add 50uL of the 1% hydrogen peroxide solution to all test wells and incubate for 10 minutes at room temperature. Do not add 1% hydrogen peroxide to standards and controls.
7. Wash the plate three times with wash buffer, invert and tap gently to remove residual buffer.
8. Dilute 50uL of the anti-Onchocerca antibody (Yellow cap) to 6 mL of antibody diluent (Blue solution). Mix thoroughly. Add 50 uL of diluted rabbit anti-onchocerca antibody to all wells. Place the plate in a humidity chamber and incubate for 45 minutes at 37°C.
9. Wash the plate three times as before.
10. Dilute 50 uL of Anti-rabbit HRPO conjugate (Purple cap) to 6 mL of antibody diluent. Add 50 uL of diluted conjugate to all wells. Place in the plate in a humidity chamber and incubate for 45 minutes at 37°C
11. Wash the plate three times as before.
12. Prepare the substrate by diluting 300uL of TBM substrate concentrate to 5700uL of Substrate Buffer. Mix thoroughly. Add 50uL to each well and incubate for 15 minutes in the dark at room temperature.
13. Stop the reaction by adding 50uL of Stopping Solution to each well (blue colour will turn yellow).
14. Read the plate using a plate reader at 450 nm, or dual wavelengths of 450/620 nm.

Figure 1: ELISA plate layout



Figure 2: Typical Standard Curve for quantification.

READING AND INTERPRETATION OF RESULTS

Titre Group Interpretation by Moore and Copeman (1990)

Read the microwell plate at dual wavelength in a compatible ELISA plate reader, blanked against air. For the test results to be accepted, standard 1 (Negative control) and Standard 7 (Positive Control) must read as follows:

- Standard 1 OD <0.2
- Standard 7 OD >2.0

A typical standard curve for the seven standard antigens is shown below. These values are used to allocate the test samples into one of eight titre groups according to the table shown below.

Using the seven standard samples it is possible to allocate the test samples into eight titre groups according to the following table. The titre groups are very useful for population studies. Moore and Copeman (1990) allocated antigen units to the seven controls.

Table 1

Allocation of samples to titre groups				
Titre Group	Absorbance	Standard No.	Antigen units	Average OD values
1	< Control sample No. 1	1	<10	0.059
2	< Control sample No. 2	2	32	0.335
3	< Control sample No. 3	3	128	0.556
4	< Control sample No. 4	4	512	0.992
5	< Control sample No. 5	5	2,048	1.520
6	< Control sample No. 6	6	8,192	2.243
7	< Control sample No. 7	7	32,000	2.809
8	> Control sample No. 7			>3.0

Note: Standards 1-7 results shall be within range of a linear Standard Curve. Refer to the Average OD values for each standard.

Standard 1 control sample does not contain any antigen (<OD 0.2). Test samples allocated in titre group 1 can be considered non-reactors (negative). Samples allocated in titre group 2 are equivocal reactors (grey area) and the test may be repeated. Test samples allocated in titre groups 3-7 are positive reactors containing circulating filarial antigen (CFA).

A previous study in New Guinea, Ghana, the Philippines and India this group represented up to 10% of the test samples. None of the Australian samples from an uninfected population were allocated to group 3. It is very likely that samples allocated to group 3 are reacting in the assay. Further data will be collected on this group to determine their status.

Recommended Low Prevalence Interpretation Method

A low prevalence interpretation is introduced to this assay due to the significant impacts of mass drug administration (MDA) aimed at eliminating filariasis. This method applies to low prevalence studies where antigen levels had significantly decreased, close to the end-point, after repeated MDA where antigen levels are below Titre Group 3.

Use Standard 2 as the reference for determining a cut-off value to define a low positive sample and a negative sample. If Standard 2 OD reading is within +/- 10% of OD 0.35, and the curve is linear or within range of the average values on Table 1, this value shall be used as the cut-off value for result interpretation. If Standard 2 is not within range of OD 0.35, a default cut-off value of OD 0.35 shall be used.

The default cut-off value of OD 0.35 was determined using 300 non-endemic serum samples from Australia. The cut-off value was determined by calculating the mean OD +3 Standard Deviation. Standard 2 is expected to give a value of OD 0.35 at 450/620nm wavelength but may vary between assays.

Studies of serum samples of uninfected Australian residents

Serum:	Samples	Mean	Standard Deviation
Study 1	303	OD 0.149	0.019
Study 2	66	OD 0.063	0.129

Dried Blood Sample on filter paper :

Study 1	100	OD 0.149	0.045
---------	-----	----------	-------

NOTES ON THE ASSAY

This kit has been designed to be used in conjunction with serum or plasma samples. Samples are boiled in an EDTA solution and centrifuged. This treatment has been shown to dissociate antigen/antibody complexes. The target antigen is heat stable and is retained in the supernatant. Following treatment there can be up to a fourfold rise in antigen titre. However, the boiling treatment can be optional, a study by Reeve, D and Melrose W, 2014 and in-house studies (unpublished) found that the boiling treatment may not significantly increase the sensitivity of the assay however, the use of an appropriate dilution may nullify the difference if samples are not treated by boiling.

For filter paper blood samples, a 1% hydrogen peroxide (H₂O₂) solution step is required to reduce the interference of endogenous peroxidase released by the haemolysis of red blood cells on the filter paper. Test samples can be handled in a 96 well microtitre format using racked tubes. The racks are modified to allow them to be placed into a boiling water bath. Subsequently the racks can be centrifuged in a rotor with microtitre buckets at approximately 2,000 g for 15 minutes.

The 96 well microtitre plates supplied are coated with a monoclonal antibody (Og4C3) which has been shown to specifically recognise only *Wuchereria bancrofti* antigen in human sera. Og4C3 will not cross react with human sera infected with *Onchocerca volvulus*, *Brugia malayi*, *Brugia timori*, *Loa loa*, *Mansonella perstans*, *Strongyloides stercoralis*, *Draconculus medinensis* or *Ascaris lumbricoides*. The cattle parasite. *Onchocerca gibsoni* is recognised by this monoclonal antibody and is used to standardise the ELISA.

The indicator antibody is produced by vaccinating rabbits with purified *Onchocerca gibsoni* antigen and the conjugate is a goat anti-rabbit globulin conjugated to horseradish peroxidase. The addition of TMB substrate solution produces a blue colour during development of positive samples, a reaction which is stopped by the addition of Stopping solution in which changes the blue colour to yellow.

A humid chamber for incubating plates can be prepared by laying a few layers of wet laboratory tissue flat at the bottom of a plastic box with a lid. Plates can be placed on top of the wet laboratory tissue and incubated in a 37°C incubator or water bath.

FILARIASIS ANTIGEN Og4C3 ELISA

ELISA Kit for the detection and quantification of *Wuchereria bancrofti* Antigen

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
References

Referencias

- More, S.J., and Copeman, D.B. (1990) A highly specific and sensitive monoclonal antibody-based ELISA for the detection of circulating antigen in bancroftian filariasis. *Tropical Medicine and Parasitology* 41: 403- 406
- Partono, F., Oernijati, S. and Hudojo, I. (1977) Malayan filariasis in Central Sulawesi (Celebes), Indonesia. *Southeast Asian Journal of Tropical Medicine and Public Health* 8: 452-458
- Reeve, D. and Melrose, W. (2014), Evaluation of the Og4C3 Filter Paper Technique in Lymphatic Filariasis Prevalence Studies, *Lymphology* 47 (2014) 65-72.
- Turner P., Copeman B., Gerisi, D. and Speare R (1993) A comparison of the OG4C3 antigen capture ELISA, the Knott test, an IgG4 assay and clinical signs, in the diagnosis of Bancroftian filariasis. *Tropical Medicine and Parasitology* 44: 45-48
- World Health Organisation (1984) Lymphatic Filariasis. Fourth Report of the WHO Expert Committee on Filariasis, Geneva. WHO Technical Services 702

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