MALARIA Ag CELISA

English
Product Code:
KM2/KM2BP

INTENDED USE

The Malaria Ag CELISA™ has been designed as a confirmatory test for *Plasmodium falciparum* (P.f) malaria in situations where traditional diagnosis is unclear, for screening blood transfusion products, or to confirm cases of travel-related infection. It is **not** intended to replace the conventional blood film diagnosis.

INTRODUCTION

Malaria is still one of the most prevalent diseases in the world with a major impact on developing tropical countries. Over 2 billion people live in malarious areas of the world and over 200 million are infected each year. Malaria is caused by protozoan parasites of the genus *Plasmodium*. Four species of the parasite commonly infect man, the most important being *Plasmodium falciparum* which is responsible for approximately 3 million deaths per year.

The Malaria Ag CELISA™ kit is a monoclonal antibody-based assay specific for P.f malaria. The assay detects HRP-2, an antigen produced at the merozoite (red blood cell) stage of the infection. This antigen circulates in the blood for up to 14 days post infection.

Malaria should be considered in the diagnosis of **any** fever of unknown aetiology in persons who have travelled to malarious areas, have recently received whole-blood transfusions or who live or work near airports. There are currently 91 countries and territories with malaria transmission. Acute malaria may be mistaken for tuberculosis, yellow fever, appendicitis, influenza, urinary tract infection or other diseases.

PRINCIPLE OF THE TEST

The Malaria Ag CELISA™ is a sandwich ELISA where microwells used in the assay are pre-coated with anti-P.f monoclonal antibody (MAb). The addition of a sample in the well allows the binding of anti-P.f MAb and HRP-2 complex if the protein is present in the sample. A wash step removes the unbound components. A second MAb conjugated to enzyme horseradish peroxidase is added as the detector which binds to the complex. The addition of a substrate solution develops the bound complex to a blue colour. Stopping solution is added to terminate the reaction, turning the microwells from blue to a yellow colour in which intensity is proportional to the HRP-2 in the sample.

CONTENTS OF THE KIT

The Malaria Ag CELISA is available in 2 different formats:

		KM2	KM2BP
MAMW	CELISA Plate – 1x96 wells - (single use only)	2 plates	10 plates
CONTROL1	Positive Control (freeze dried)	1 x 0.5mL	2 x 0.5mL
CONTROL2	Negative Control	1 x 2.5mL	1 x 5.0mL
MAPO	Enzyme Conjugate (200x)	1 x 0.12mL	1 x 0.7mL
MACD	Conjugate Diluent	1 x 24mL	1 x 120mL
MAPT	PBS/Tween (20x)	1 x 125mL	2 x 250mL
MASC	Substrate Chromogen (TMB) (20x)	1 x 1.2mL	1 x 7mL
MASB	Substrate Buffer	1 x 24mL	1 x 125mL
MASS	Stopping Solution	1 x 12mL	2 x 30mL

Store all components at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Malaria positive blood samples, distilled water, micropipettes and tips, clean glassware or plastic containers for solutions, humid chamber, ELISA washer, spectrophotometer to read absorbances at a single wavelength of 450nm, or at dual wavelengths of 450nm and 620nm.

PRECAUTIONS

For in vitro diagnostic use only. Reagents should not be used after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Thimerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eyes and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to the Safety Data Sheet.

INSTRUCTIONS FOR USE

Preparation of Wash Buffer

If crystals are present, warm the concentrate to dissolve. For each microplate, add 50mL PBS-Tween concentrate **MAPT** to 950mL of distilled water. Label the bottle **WASH BUFFER**. Store at 2-8°C.

Preparation of Samples

Collect patient blood by standard venepuncture procedure using an anticoagulant. Lyse blood by freezing, use the lysed blood as the test specimen (**SAMPLE**). Serum or plasma may be used as an alternative to whole blood but the use of these samples may result in the loss of sensitivity. Blood samples should be stored below -10°C if the analysis is delayed.

Preparation of Positive and Negative Controls

The Negative Control **CONTROL1** is RPMI, add 100µL directly into the well as a negative control, RPMI can also be used to dilute the Positive Control. The Positive Control **CONTROL2** is a culture supernatant containing HRP-2, a known positive blood can also be used as an internal positive control. **IMPORTANT: The kit Positive Control is freeze dried at a concentration of 100 ng/mL. Reconstitute by adding 500µL of distilled water and standing for 5 minutes before mixing thoroughly. It is recommended to aliquot reconstituted Positive Control into smaller volumes to avoid repeat freeze-thaw cycles. Once reconstituted, store vials at -20°C until use. Dilute further to 1:10 in a microtube before using as a control in the ELISA.**

Assay Procedure

- Bring all reagents to room temperature (18-25°C) before use.
- Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
- Remove required number of **MAMW** strips. Reseal the foil bag containing unused microwell strips immediately.
- Reconstitute the Positive Control **CONTROL1** using the recommended method described in Page 1, **Preparation of Positive and Negative Controls**. Pipette 100µL of the **SAMPLE**, Negative Control **CONTROL2** and Positive Control **CONTROL1** (reconstituted and diluted to 1:10) into individual microwells. Test controls in duplicate. Cover and incubate for one (1) hour at room temperature in a humid chamber.

The set of instructions below is for one 8-well strip only. Calculate the required reagent for your test, allowing 1mL per 8-well strip.

- In the last 10 minutes of the incubation period, prepare the working strength **CONJUGATE**. For each 8-well strip, add 5µL of conjugate concentrate **MAPO** to 995µL of conjugate diluent **MACD** and mix thoroughly.
- Wash the wells 4 times preferably using an automatic plate/strip washer or manually as follows:
 - Empty contents from the wells. Refill with the **WASH BUFFER**.
 - Repeat this process a further four (4) times. After the fifth wash, tap dry the plate on absorbent tissue.
 - NB: take care when flicking out plates, hold side of frame firmly to hold strips in place.
- Add 100µL of **CONJUGATE** to each well. Incubate for one (1) hour at room temperature in a humid chamber.
- In the last 10 minutes of the incubation period, prepare the working strength **SUBSTRATE**. For each 8-well strip, add 50µL of substrate chromogen **MASC** to 950µL of substrate buffer **MASB** and mix thoroughly. The stability of the solution is 30 minutes.
- Repeat washing as in step 6.
- Add 100µL of fresh **SUBSTRATE** and incubate in the dark (covered) at room temperature for 15 minutes.
- Add 50µL of Stopping Solution **MASS**. Tap the plate to mix.
- Read the results visually or in a spectrophotometer at 450nm, or 450nm/620nm, blanking the machine on air.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

Visually

Observe the colour intensity of the control and specimen wells. The Positive Control should be blue before, and yellow after stopping.

Photometrically

Read the microwell plate at 450nm or 450nm / 620nm in a compatible ELISA plate reader, blanked against air.

For the test results to be accepted the Negative Control must read as follows:

	O.D Value (450nm, 450/620nm)
Negative Control	< 0.1
Positive Control	> 1.5
Cut-Off level	= Negative Control OD + 0.1

Negative blood samples should give optical density readings below 0.1 OD units. However, to allow for inter-laboratory variation we strongly recommend that each laboratory run a number of known negative blood samples to allow standardisation of the CELISA positive / negative cut-off level.

All specimens with an absorbance value above the cut-off level should be considered positive for *P. falciparum* antigen. A positive result indicates the presence of *P. falciparum* HRP-2 antigen. This is suggestive of current or very recent infection. The assay has been shown to detect *P. falciparum* infection at parasitaemia as low as 0.001%. The intensity of colour is not proportional to the level of parasitaemia. *P. vivax*, *P. ovale* and *P. malariae* infections are not detected. Please note that the test may remain positive for several days after parasites are no longer detectable in blood films.

Quantification of HRP-2

An improved version of the Malaria Ag CELISA kit was developed specifically for quantification of Pf-HRP-2, the Quantimal™ CELISA Ultra-sensitive Pf-HRP-2 Malaria (Product Code: KM8). The estimated lower detection limit of KM8 is as low as 40 pg/mL. Please contact your distributor or Cellabs for more information on the suitable ELISA kit for your quantification needs.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE MALARIA AG CELISA

Refer to summary table at end of insert. All data on the Malaria Ag CELISA can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying causative agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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Cellabs Pty Ltd
WMDE B.V