Quantimal™ pLDH Malaria CELISA

English
Product Code: KM7/KM7BP

BACKGROUND AND INTENDED USE

The Quantimal™ pLDH Malaria CELISA is for the detection of malaria parasite lactate dehydrogenase (pLDH), a biomolecule common to all *Plasmodium spp.*, used as a marker for malaria diagnosis. It is well established that pLDH in malarial blood samples is indicative of active infection and its reduction mirrors the clearance of infection after treatment. Unlike histidine rich proteins like HRP-2, the pLDH enzyme is secreted only by live parasites and survives briefly in the patient blood. This kit can be used for monitoring the resistance of *Plasmodium* parasites to anti-malarial chemotherapy. It can also be used as a quantitative assay, for screening malaria in blood banks and routine clinics, and for in-vitro evaluation of anti-malarial drugs. This kit is a confirmatory test, not intended to replace the conventional blood film diagnosis, to be used with a variety of clinical samples such as whole blood, patient blood culture, plasma, serum or in-vitro preparations for drug resistance studies.

The sandwich ELISA principle is employed using microwells pre-coated with anti-pLDH monoclonal capture antibody. A second anti-pLDH monoclonal, the detector antibody, is labelled with horseradish peroxidase enzyme. If pLDH is present in a sample, it binds to the antibody coated on the well. A washing step removes the unbound components to allow a complex to form once the labelled detector antibody is added in an incubation step. A blue colour complex is developed with the addition of a substrate TMB which turns yellow after the addition of a stopping solution. The colour intensity is proportional to the amount of pLDH detected in a sample.

CONTENTS OF THE KIT

		KM7	KM7BP
MPMW	Celisa plate – 2 x 96 wells (single use)	2 plate	10 plates
CONTROL +	Positive Control	1 x 0.5 mL	1 x 1.0 mL
CONTROL -	Negative Control	1 x 2.5 mL	1 x 5.0 mL
MPPQ	Enzyme Conjugate (200x)	1 x 0.12 mL	1 x 0.7 mL
MPCD	Conjugate Diluent	1 x 24 mL	1 x 120 mL
MPPT	PBS/Tween (20x)	1 x 125 mL	2 x 250 mL
MPSC	Substrate Chromogen (TMB) (20x)	1 x 1.2 mL	1 x 7.0 mL
MPSB	Substrate Buffer	1 x 24 mL	1 x 125 mL
MPSS	Stopping Solution	1 x 12 mL	2 x 30 mL

Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Distilled water, micropipettes and tips, clean glassware or plastic containers for solutions, humid chamber, ELISA washer, ELISA plate reader to read absorbances at a single wavelength of 450nm or at dual wavelength of 450nm and 620nm

PRECAUTIONS

For in-vitro diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. Thiomerosal preservative added to some components is a poison. Exercise caution when handling these components. The stopping solution is corrosive. Avoid contact with skin, eye and mucous membranes. Dispense all reagents with care to avoid cross contamination of wells. Avoid exposure of the substrate to light. Treat all clinical and control material as though potentially infectious and dispose of in accordance with local operating regulations. For further information, please refer to Material Safety Data Sheet.

INSTRUCTIONS FOR USE

IMPORTANT: The positive control is time and temperature sensitive and stable for a maximum of 4 months at 2-8°C. We strongly recommend that the vial is stored frozen at -20°C to -80°C immediately when kits are received. The vial can remain frozen until use.

Preparation of Wash Buffer

Wash buffer is provided at 20x concentrate. If crystals are present, pre-warm the concentrate and stir on a magnetic stirrer until all crystals are dissolved. For each microplate, add 50mL PBS-Tween concentrate **MPPT** to 950mL of distilled water. Label the bottle 1x **WASH BUFFER** and store at 2-8°C. This buffer is used for diluting test samples, diluting the conjugate concentrate **MPPQ** and washing the microplates.

Preparation of samples

Collect patient blood by standard venepuncture procedure using an anticoagulant. Blood specimens should be stored below -10°C if the analysis is delayed. Fresh, refrigerated or frozen samples of serum, plasma or blood may be used. Avoid contamination by collecting aseptically. The **CONTROL +** is culture supernatant and needs to be diluted 1:1 with the **CONTROL -**. This can be prepared directly into the well by adding 50µL of **CONTROL +** + 50µL of **CONTROL -**, a total volume of 100µL each well. The negative control **CONTROL -** is RPMI and this is to be tested neat, directly into the well. All test samples are to be diluted 1:1 with prepared diluent (1x **WASH BUFFER**) and can be prepared directly into the wells.

Assay Procedure

- Bring all reagents to room temperature (18-25 °C) before use.
- Prepare **WASH BUFFER** (see Preparation of Wash Buffer)
- Remove required number of **MPMW** strips. Firmly reseal the snap-lock foil bag containing unused microwell strips.
- Controls: Pipette 100µL volume duplicates of each of the controls according to the above recommended instructions.
- Samples: Pipette 50µL of 1x **WASH BUFFER** directly into the well(s) followed by 50 µL of each of the sample(s), taking care not to cross-contaminate the neighbouring wells. Change pipette tips after each sample addition. Gently tap the side of the ELISA frame to give the wells a gentle mix.
- Place the plate in a humidity chamber and incubate for 30 minutes at 37°C. A humidity chamber can be any plastic box with a flat bottom and a lid.
- In the last 10 minutes of the incubation period, prepare the working strength **CONJUGATE**. Add 5µL of conjugate concentrate **MPPQ** to 995µL of conjugate diluent **MPCD** and mix thoroughly (allow 1mL per strip of 8 wells).
- Wash the wells four times, preferably using an automatic plate/strip washer or manually as follows:
 - Empty contents from the wells. Refill with the **WASH BUFFER**.
 - Repeat this process four (4) times. After the fourth wash, bang inverted wells dry on absorbent tissue.
 - NB: take care when flicking out plates, hold side of frame firmly to hold strips in place.
- Add 100µL of diluted **CONJUGATE** to each well. Incubate for 30 minutes in a humid chamber at 37°C.
- In the last 10 minutes of the incubation period, prepare the working strength **SUBSTRATE**. Add 50µL of substrate chromogen **MPSC** to 950µL of substrate buffer **MPSB** and mix thoroughly (allow 1mL per strip of 8 wells). The stability of the solution is 15 minutes.
- Repeat washing as in step 8.
- Add 100µL of fresh **SUBSTRATE** and incubate in the dark (covered) at room temperature for 15 minutes.
- Add 50µL of Stopping Solution **MPSS**. Tap the plate to mix.
- Read the results visually or in a spectrophotometer preferably at 450nm/620nm (or 450nm).

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

Visually

Observe the colour intensity of the control and specimen wells. The Positive Control should be blue before and yellow after stopping.

Photometrically

Read the microwell plate at 450nm or 450/620 in a compatible ELISA plate reader, blanked against air. For the test results to be accepted the controls must read as follows:

	O.D. Value (450nm or 450/620nm)
Negative Control	O.D < 0.2
Positive Control	O.D > 1.5
Cut-off value (COV)	O.D 0.25

The cut-off value (COV) specified for the Quantimal™ pLDH Malaria CELISA was determined by testing 292 samples from a non-endemic area. Any sample with an optical density (OD) value above the COV is considered positive for pLDH. A positive result indicates the presence of pLDH from *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* or *Plasmodium knowlesi*. The Quantimal™ pLDH Malaria CELISA does not distinguish between different species of *Plasmodium* infections. To detect specific infection caused by *Plasmodium falciparum*, contact Cellabs Pty Ltd for information on the Malaria Antigen (HRP-2) CELISA.

Quantification

This assay can be used for the quantification of pLDH, please contact Cellabs at enquiries@cellabs.com.au for further details on how to construct a standard curve using the positive control provided with the kit.

WASTE DISPOSAL

Dispose of any unused component as bio-hazardous waste. For more information, please refer to the MSDS.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims and performance of the kit. Cellabs and its agents or distributors are not liable for damages under these circumstances.

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WMDE B.V.