FIGURE 1: PNEUMO CEL INDIRECT DIAGRAM FOR USE


TABLE 1: SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE PNEUMO INDIRECT CEL

Trial	Sensitivity	Specificity	Repeatability	Reproducibility
A	100%	100%	-	-
B	100%	100%	-	-
C	100%	100%	-	-
D	100%	95.1%	-	-
E	-	-	100% Correlation	100% Correlation
<b>Not cross reactive with</b> <i>Cryptococcus neoformans</i> , <i>Candida tropicalis</i> , <i>Saccharomyces cerevisiae</i> , <i>Aspergillus niger</i> , <i>Candida krusei</i> , <i>Candida albicans</i>				

EXPLANATION OF SYMBOLS

- Consult Instructions for Use
- In Vitro Diagnostic Medical Device
- Temperature Limitation  
2°C - 8°C
- Batch
- Control Positive
- Use By/Expiration Date
- Do Not Re-use

Cellabs Pty Ltd  
Unit 7, 27 Dale Street  
Brookvale, NSW 2100 Australia  
Tel: +61 2 9905 0133 Fax: +61 2 9905 6426  
Web: <http://www.cellabs.com.au>  
Email: [sales@cellabs.com.au](mailto:sales@cellabs.com.au)

WMDE  
Bergerweg 18  
6085 AT Horn  
The Netherlands

LP2.4.1  
23 May 2017

PNEUMO CEL  
INDIRECT IF TEST

Product Code: KP2

INTENDED USE

The **Pneumo Cel Indirect IF Test** is an *in vitro* indirect Immunofluorescence test for the detection of *Pneumocystis carinii* (*P.jirovecii*) in human specimens.

INTRODUCTION

*P. carinii* (*P.jirovecii*) (PC) is an extracellular eukaryotic organism shown to be related to the fungi (1). It is an opportunistic pathogen causing interstitial pneumonia in immunosuppressed hosts, such as premature infants, children with congenital immunodeficiencies and organ transplant recipients (2,3,4). With the advent of the acquired immunodeficiency syndrome (AIDS), *P. carinii* (*P.jirovecii*) pneumonia (PCP) has emerged as the most common life-threatening opportunistic infection in these patients. PCP occurs in over 75% of AIDS patients (5,6,7) and is usually associated with diminished CD4+ cells (8). PC is generally confined to the lungs, although there have been recent reports of disseminated extrapulmonary disease (9) which may be as a result of the common use of aerosolised pentamidine (10).

The most common form of PC is the thick-walled cyst (5-8µm) containing up to 8 intracyclic bodies (sporozoites, 1-2µm). These mature into extracyclic pleomorphic forms, the trophozoites (2-5µm). The development cycle occurs in both the alveolar epithelial cells and in the alveolar spaces.

PRINCIPLE OF THE TEST

PC may be identified by a variety of histochemical stains such as methenamine silver, Giemsa and toluidine blue (11). The stains can, however, be difficult to perform, they are non-specific and they may not detect all stages of the life cycle of the organism. Experience is also required for interpretation. More recently, immunofluorescence assays (IFA) using monoclonal antibodies to PC have been developed. These techniques are rapid, easy to perform and specimens from induced sputum and bronchoalveolar lavage can be used, thus avoiding the need for more invasive procedures such as lung and transbronchial biopsies. IFA has been shown to be both specific and sensitive (12,13), and also more sensitive than histochemical staining (14).

The **PNEUMO-CEL INDIRECT IF TEST** is an indirect immuno-fluorescence test for PC organisms in methanol-fixed specimens. The fluorescein-labelled mouse monoclonal antibodies bind specifically to the murine raised monoclonal, initially bound to the organisms. Washing steps, following both the monoclonal and fluorescein conjugated antibodies, removes any unbound antibody. When viewed under a fluorescence microscope, PC are seen as bright, apple green bodies which may be present singularly or in cluster. Trophozoites, cysts and the extra cellular matrix, which binds cysts, are stained.

CONTENTS OF THE KIT

		KP2 Standard	Bulk
	Monoclonal Antibody Reagent	1.25mL	5mL
	Anti Mouse FITC Reagent	1.25mL	5mL
	Mounting Fluid	2.5mL	2 x 2.5mL
	Tests	50	200

Materials are supplied ready for use. Store at 2-8°C. Expiry dates are clearly marked on each kit component and on the box. Expiry dates do not change once opened.

MATERIALS REQUIRED BUT NOT PROVIDED

Positive control human PC; microscope slides with 6-8mm diameter wells; precision pipette for delivering 25µL; methanol for specimen fixation; humid chamber; wash bath; phosphate buffered saline (PBS) for washing step; cover slips; non-fluorescing immersion oil; and fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x400-x1000 magnification.

PRECAUTIONS

For *in vitro* diagnostic use only. Do not use after the expiry date shown on the label. If protective packaging is damaged, contact your local distributor and ask for a replacement. Do not mix reagents from different kits. The Pneumo Cel Reagent has been optimised for use with Cellabs Mounting Fluid. Evans Blue dye contained in the test reagent is a possible carcinogen, therefore avoid contact with the skin. Patient specimens and the positive control slide should be handled as though potentially infectious. A Positive Control Slide must be run with each test run. Refer to Material Safety Data Sheet (MSDS) for further information.

SPECIMEN COLLECTION

Clinical specimens should be processed as soon as possible after collection. Induced sputum, bronchial wash, bronchoalveolar lavage and tissue impression smears may be used. The presence of mucous in specimen may prevent adequate staining. Work in a biohazard hood. A positive control slide should be included with each series of test specimens.

Preparation of Induced Sputum Specimens

Prepare 0.3% dithiothreitol (DTT) in distilled water. Mix equal volumes of sputum and DTT solution. Mix vigorously (vortex) to dispense clumps. Incubate for 30 minutes at 37°C. Centrifuge for 5 minutes at 1500g. Decant supernatant. Resuspend pellet in 0.2-0.5mL of PBS and mix vigorously.

Preparation of Bronchial Wash and Bronchoalveolar Lavage Specimens

Centrifuge the specimen at 1500g for 10 minutes. Remove the supernatant leaving 0.5mL. Decant supernatant. Resuspend pellet in 0.2-0.5mL of PBS and mix vigorously

PREPARATION OF SLIDES

Apply 25µL to 1 well (0.8cm) slide and air dry. Fix slide for 5 minutes in acetone. If slides are not stained immediately they may be stored below 0°C for several weeks.

INSTRUCTIONS FOR USE

1. Add 25µL of monoclonal reagent to the fixed test and control slide, covering well area.
2. Incubate the slides in the dark at 37°C in a humid chamber for 30 minutes. Do not allow slides to dry out.
3. Wash the slides gently in PBS bath, drain the slides and remove excess moisture until just dry.
4. Add 25µL of the anti-mouse Ig FITC reagent to the specimen. Incubate slides in the dark as in step 2.
5. Wash the slides gently in PBS bath, drain the slides and remove excess moisture until just dry
6. Add a drop of to the slide well. Place a coverslip on top of the drop and remove air bubbles.
7. Scan the entire specimen using a fluorescence microscope initially at X400 magnification, then at X1000 for confirmation. Read immediately or store at 2-8°C in the dark for up to 24 hours.

READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

1. PC should exhibit bright, apple-green fluorescence of the round, thick-walled cysts and pleomorphic trophozoites in contrast to the reddish-brown colour of counterstained material. The organisms occur either singularly or more commonly in clusters in an extracellular matrix, which is also stained. A test is considered positive if two or more cysts are present.
3. Negative specimens should not display fluorescent organisms of characteristic morphology.
4. Fluorescent material which is not apple green and can be distinguished from typical PC forms should be disregarded.

WASTE DISPOSAL

Dispose of any unused components as biohazardous waste. Where the test reagent has been disposed of in the sink, ensure it is flushed with large quantities of water (as the sodium azide it contains may react with copper/lead plumbing systems). For more information, please refer to the MSDS.

SENSITIVITY, SPECIFICITY, & OTHER DATA ON THE PNEUMO CEL

Refer to summary table at end of insert. All data on the Pneumo Cel can be obtained in the product information sheet. Please ask your local distributor or contact Cellabs.

INDEMNITY NOTICE

Modifications or changes made in the recommended procedure may affect the stated or implied claims. A positive or negative result does not preclude the presence of other underlying respiratory agents. Cellabs and its agents and distributors shall not be liable for damages under these circumstances.

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



Images of human *Pneumocystis carinii* stained with Pneumo Cel Indirect IFA.

### REFERENCES

- Edman, J.C., Kovacs, J.A., Masur, H., Santi, D.V., Elwood, H.J., Sogin, M.L. (1988). Ribosomal RNA sequence shows *Pneumocystis carinii* to be a member of the fungi. *Nature*, 334:519-22.
- Walzer, P.D., Schultz, M.G., Western, K.A., Robbins, J.B. (1976). *Pneumocystis carinii* pneumonia and primary immune deficiency diseases. *Nat. Cancer Inst. Monogr.* 43:65-72.
- Gryzan, S., Paradis, I.L., Zeevi, A., Duquesnoy, J., Dummer, S., Griffith, B.P., Hardesty, R.L., Trento, A., Nalesnik, M.A., Dauber, J.H., (1988). Unexpectedly high incidence of *Pneumocystis carinii* infection after lung heart transplantation. *Am. Rev. Respir. Dis.* 137:1268-1274.
- Hardy, A.M., Wajszczuk, C.P., Suffrendini, A.F., Hakala, T.R., Ho, M., (1984). *Pneumocystis carinii* Pneumonia in renal transplant recipients treated with cyclosporine and steroids. *J. Infect. Dis.* 149:143-147.
- Hughes, W.T., (1987) *Pneumocystis carinii* pneumonia. *N. Eng. J. Med.* 317:1021-3.
- Kovacs, J.A., Masur, H. (1988). *Pneumocystis carinii* pneumonia: therapy and prophylaxis. *J. Infect. Dis.* 158:254-9.
- Murray, J.F., Garay, S.M., Hopewell, P.C., Mills, J., Snider, G.L., Stover, D.E. (1987). Pulmonary complications of the acquired immunodeficiency syndrome: an update: report of the Second National Heart, Lung and Blood Institute Workshop. *Am. Rev. Respir. Dis.* 135:504-509.
- Phair, J., Munoz, A., Detels, R., Kalslow, R., Rinaldo, C., Saahr, A. (1990). The risk of *Pneumocystis carinii* pneumonia among men infected with the human immunodeficiency virus type 1. *New Engl. J. Med.* 322:161-165.
- Pilon, V., Echols, R., Celis, J., Elmendorf, S., (1987). Disseminated *Pneumocystis carinii* pneumonia infection in AIDS. *N. Eng. J. Med.* 316: 1410-1411.
- Raviglione, M.C. (1990). Extrapulmonary pneumocystosis: the first 50 cases. *Rev. Infect. Dis.* 12:1127-1138.
- Ruskin, J. (1996). Newer developments in diagnosis and treatment of *Pneumocystis* infections. *Curr. Clin. Top. Infect. Dis.* 7:194-215.
- Kovacs, J.A., Ng, V.L., Masur, H., Leoung, G., Hadley, W.K., Evans, G., Lane, H.C., Ognibene, F.P., Shelhamer, J., Parillo, J.E., Gill V.J. (1988). Diagnosis of *P. Carinii* pneumonia: improved detection in sputum with use of monoclonal antibodies. *New Eng. J. Med.* 318:589-593.
- Wolfson, J.S., Waldron, M.A., Sierra, L.S. (1990). Blinded comparison of a direct immunofluorescent monoclonal antibody staining method and a Giemsa staining method for identification of *P. Carinii* in induced sputum and bronchoalveolar lavage specimens of patients infected with human immunodeficiency virus. *J. Clin. Micro.* 28:2136-2138.
- Ng, V.L., Virani, N.A., Chaisson, R.E., Yajko, C.D., Sphar, H.T., Cabrian, K.C., Rollins, N., Charache, P., Krieger, M., Haldey, W.K., Hopewell, P.C. (1990). Rapid detection of *Pneumocystis carinii* using a direct fluorescent monoclonal antibody stain. *J. Clin. Micro.* 28:2228-2233.