


Leishmania CEL IFA

The Cellabs Leishmania CEL IFA kit is a direct immunofluorescence assay designed to detect all Leishmania protozoa: promastigotes from sand-flies or culture and amastigotes from human biopsy material.

Convenient & easy to use

Single 30 minutes incubation
Ready to use reagents (50 tests)
Control slide included in kit
Long shelf life (18+ months)



Leishmania promastigotes

Highly reliable clinically

Excellent fluorescence
Direct demonstration of the parasitemia

Reliable diagnostic confirmation

- Direct demonstration of the parasite
- Single reagent to detect ALL parasite stages
- Simple and rapid diagnosis in less than 1 hour
- Convenient packaging and long shelf life



Cellabs Pty Ltd

P.O. Box 421 Brookvale

NSW 2100 Australia

Ph: + 61 2 9905 0133

Fax: + 61 2 9905 6426

Website: www.cellabs.com.au

Email: sales@cellabs.com.au