

Visceral Leishmania IgG CELISA

The Cellabs Visceral Leishmania IgG CELISA is a direct ELISA assay using a novel antigen from cultivating *L. donovani* promastigotes in protein-free media. It detects parasite-specific IgG in serum.


Convenient & easy to use

Results in about 3 hours
Ready to use reagents (96 tests)
Can be read visually



Sensitive & specific

98% specificity
97% sensitivity



A clinically reliable test

- A reliable indicator for epidemiological studies
- A practical tool for blood bank screening
- Can be used for post-treatment evaluation



Cellabs Pty Ltd

P.O. Box 421 Brookvale
NSW 2100 Australia
Ph: + 61 2 9905 0133
Fax: + 61 2 9905 6426
Website: www.cellabs.com.au
Email: sales@cellabs.com.au