

Malaria Antigen CELISA

Malaria Antigen CELISA is a direct ELISA kit for the specific detection in whole blood or serum of HRP-2 secreted by *P. falciparum* at the merozoite (red blood cell) stage of the infection.

Easy to use

Easily automated protocol
Rapid results (2 hours)
Can be used in combination
with Malaria Antibody CELISA



Superior clinical reliability

More sensitive (5-50 parasites/ μ l) in practice than usually performed thick or thin film microscopy

Reliable alternative to microscopy

- Alternative to thick or thin smear microscopy
- High throughput possible for screening
- Specific for potentially lethal *P. falciparum*
- Highly sensitive (5-50 parasites/ μ l)
- 192 tests kit format with long shelf life



Cellabs Pty Ltd

P.O. Box 421 Brookvale
NSW 2100 Australia
Ph: + 61 2 9905 0133
Fax: + 61 2 9905 6426
Website: www.cellabs.com.au
Email: sales@cellabs.com.au