


Pan Malaria IgG CELISA

Pan Malaria Antibody CELISA is a direct ELISA kit for the specific detection in serum of IgG antibodies produced by *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae*.

Easy to use in blood banks

Easily automated protocol
Rapid results (2 hours)
Can be used in combination
with Malaria Antigen CELISA



Superior clinical reliability

Widely used in endemic countries blood banks
Recombinant antigens ensure reproducibility of results

Reliable in blood bank screening

- Alternative to IFAT microscopy
- Several publications confirm validity
- High throughput for screening in blood banks
- Reduce blood donation rejection in Malaria endemic countries
- 192 tests kit format with long shelf life



Cellabs Pty Ltd

P.O. Box 421 Brookvale
NSW 2100 Australia
Ph: + 61 2 9905 0133
Fax: + 61 2 9905 6426
Website: www.cellabs.com.au
Email: sales@cellabs.com.au