

T. cruzi IgG CELISA

The Cellabs *T. cruzi* IgG CELISA is a direct ELISA kit designed to detect antibodies produced during *T. cruzi* infections (Chagas disease).

Blood bank screening

Easily automated protocol
Most reagents ready to use
192 tests pack size
Long shelf life



Clinically reliable

Optimised for use in Latin America
Native antigen kit providing superior sensitivity

Reliable blood bank screening

- **98% sensitivity and 98% specificity**
- **Designed for blood bank screening**
- **Easily automated protocol**
- **Visual or photometric detection**
- **Results in about 3 hours**



Cellabs Pty Ltd

P.O. Box 421 Brookvale
NSW 2100 Australia
Ph: + 61 2 9905 0133
Fax: + 61 2 9905 6426
Website: www.cellabs.com.au
Email: sales@cellabs.com.au